



Philips DreamWear
Mask

3 cushion sizes
Medium frame



HH1124

Like wearing nothing at all*

With the freedom and comfort of DreamWear Gel Pillows, you can sleep how you want.

Innovative design

- Air tube connection at the top of the head
- Gel pillows cushion
- Supports your sleep therapy
- Modular design

Comfortable design

- Easy to fit
- Soft fabric wraps
- Slim-line headgear
- Soft, flexible frame

Open field of vision

- Open design offers an extremely wide field of vision

PHILIPS

Mask
3 cushion sizes Medium frame

HH1124/00

Highlights

Easy to fit



Choose the cushion so the pillows cushion blue gel bases sit comfortably under the nose without any gaps. Only the soft silicone tip should be inside your nostrils. Do not push any part of the blue gel base inside your nose.

Sleep how you want



By moving the tube location to the top of the head, DreamWear Gel Pillows offers more freedom of movement during the night. Users stated they didn't need to choose a sleep position when wearing DreamWear Gel Pillows.

Gel pillows cushion



Less discomfort on the face, nose and nostrils. Users rated DreamWear Gel Pillows more comfortable than their prescribed mask.

Modular design



Modular design for quick assembly and cleaning, with interchangeable pillows and nasal cushions.

Slim-line headgear



DreamWear combines comfort with aesthetics. Users reported that the DreamWear Gel Pillows mask came closest to making them feel like they have nothing on their face versus their prescribed mask.

Soft fabric wraps



Provides a soft interface to the cheeks.

Soft, flexible frame



The flexible material through which the air flows is applied gently to the skin and is pleasant to wear.

Unobstructed view



Open design makes it easy to wear glasses, read, watch TV, and use a computer or tablet before falling asleep.

Supports your sleep therapy



DreamWear gives you the freedom to choose the cushion that works best for you, giving you freedom of movement and better quality of sleep.

Specifications

Product details

- Number of parts: 5
- Materials: Nylon, Polyurethane foam, Polycarbonate, Polyester, Polyurethane gel, Silicone, Spandex
- Operating pressure: 4 to 20 cmH2O
- Sound level: 18 dBA
- Warranty: 90 days

Package contents

- Include: Small Cushion, Medium Cushion, Large Cushion, Medium Frame

Cleaning: Mask

- Hand wash with mild detergent: Rinse thoroughly, air dry¹

Cleaning: Headgear, tubing

- Hand wash with mild detergent: Rinse thoroughly, hang to dry¹

Intent of use

- To be used on patients >30kg: with CPAP prescription

Replacement

- Inspect daily for wear: Replace at first sign of wear



Issue date 2018-06-18

Version: 5.0.3

12 NC: 0000 011 24984
UPC: 6 06959 02678 0

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* *2017 Philips User Preference Questionnaire

* The scope of supply includes a standard mask frame (size M) and three cushions for all sizes (S, M, L). The sleep therapy device is not included.

* The costs of DreamWear are private and are usually not reimbursed by the statutory health insurance.