




Philips DreamWear  
Full Face Cushion

Full Face Replacement Part  
Medium cushion



HH1128

## Like wearing nothing at all.\*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

### Innovative design

- A modular design for more options

### Comfort and ease of use

- Easy to use
- Under the nose full face cushion

**PHILIPS**

# Specifications

## Product details

- Materials: Polycarbonate, Silicone
- Warranty: 90 days

## Cleaning: Cushion

- Hand wash with mild detergent: Rinse thoroughly,

air dry<sup>1</sup>

## Replacement

- Inspect daily for wear: Replace at first sign of wear

# Highlights

## Easy to use



Only a few parts for easy assembly, disassembly and cleaning.

## Freedom of choice



The modular design allows users to choose the cushion that works best for them, all on one frame.

## No more nose bridge red marks



The under the nose innovative design of the cushion means no more red marks, discomfort, and irritation on the nose bridge.\*\*



Issue date 2018-02-01

Version: 1.2.1

12 NC: 0000 011 33431  
UPC: 6 06959 41458 7

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\*\* Data analysis of Aug 2017 Patient Preference trial where n = 85 and prescribed masks include Resmed Airfit F10 and F20, Respirationics Amara View and F&P Simplus; based on 10 days of use.

\* Please refer to Instructions For Use for more detailed cleaning instructions.