



Philips Nuance Pro
Mask

Gel

Pack
3 pillows (S, M, L)

A new choice for better sleep

Nuance's gel nasal pillow technology delivers exceptional comfort and seal, helping you to sleep better.

Easy to use

- Lightweight, flexible tubing for less pull on the mask.

Gel pillows

- The pillows fit different nostril sizes.

Greater stability

- Stays in place throughout the night.



HH1101

PHILIPS

Mask
Gel Pack, 3 pillows (S, M, L)

Specifications

Service

- Warranty: 90 days

Design

- Easy to use: Easy assembly and disassembly

Intent of use

- To be used on patients >30 kg (>66 lbs): with CPAP prescription

Specifications

- Operating pressure: 4–20 cmH₂O
- Sound pressure level: 19 dBA

- No. of parts: Four
- Cushion material: Silicone rubber
- Pillows/frame pads material: Urethane gel
- Headgear material: Nylon UBL/Urethane foam/Lycra
- Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry
- Cleaning: headgear, tubing: Hand wash with mild detergent
- Frame material: Thermoplastic polyester elastomer
- Swivel/Tubing/Elbow material: Polypropylene

HH1101/00

Highlights

Lightweight and simple



Easier to assemble, disassemble and use than the leading pillows mask.

Provides comfort and seal



Nuance gel pillows cushion is designed differently to provide a more comfortable experience. The gel base of the pillows should fit comfortably under your nose.

Innovative non-slip headgear



Ideal for patients who want to sleep freely.



Issue date 2018-05-23

Version: 8.1.2

12 NC: 0000 011 05168
UPC: 6 06959 02147 1

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* The purchase of this medical device does not give the right to obligatory NHS reimbursement