

Philips DreamWear
Elbow

with swivel

HH1122



Quick-release elbow with top of the head design

The quick-release elbow allows the user to easily disengage the mask from the CPAP hose without removing the mask. If sleeping positions change during the night, the swiveling connection will allow the mask to move freely.

Easy to use

- Fits the circular opening on the mask frame

PHILIPS

Elbow
with swivel

Specifications

Service

- Warranty: 90-days

HH1122/00

Highlights

Comfort and performance



Once in place, the design of the elbow directs the hose over the user's head, reducing the tug on the mask. Contains a built in exhalation port.



Issue date 2018-04-03

Version: 3.0.1

12 NC: 0000 011 16748
UPC: 6 06959 40256 0

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com