



Philips DreamStation
CPAP with Humidifier

Colorful, easy-to-use display

Display and track nightly use
Designed to reduce noise
Bluetooth enabled connectivity



HH1460/00

Take charge with ease and confidence

With easy-to-navigate menus, a sleek, compact design and quiet operation, DreamStation makes it simple to start, customize and continue your obstructive sleep apnea (OSA) therapy.

Performance

- Support for the long term
- Make progress every day
- Small and lightweight
- Performance and appeal

PHILIPS

CPAP with Humidifier

Colorful, easy-to-use display Display and track nightly use, Designed to reduce noise, Bluetooth enabled connectivity

Specifications

Product details

- Bluetooth: Built-in

Compatibility

- For use with: DreamMapper

Performance data

- Pressure mode: CPAP

- Flex mode: C-Flex
- Ramp time: 0-45 min (5-min increments)
- Starting ramp pressure: SmartRamp


Operating specifications

- Operating pressure: 4 to 20 cmH₂O
- Altitude compensation: Automatic

HH1460/00

Highlights

Performance and appeal




DreamStation offers a sleek, low-profile design. Its vibrant color screen displays daily progress and features easy-to-use menus.

Small and lightweight




DreamStation's top-mounted humidifier hose swivel and front-facing display make it easy for you to operate while lying down or sitting up in bed. And, it's small and light, making it easy to pack for travel.

Daily Progress Feedback



Each morning you are greeted with Daily Progress Feedback - a summary of how far you've come to motivate you to stick to your therapy. The intuitive, color display shows a simple trend of their nightly hours of use, followed by a summary of the last 30 days of good nights with more than four hours of use.

Bluetooth enabled



Built-in Bluetooth makes connectivity to the DreamMapper self-management tool easy. DreamMapper provides goal-setting tools, helpful videos and ongoing feedback to keep you motivated and engaged.



Issue date 2018-08-30

Version: 2.0.1

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com