

Philips InnoSpire Air Inlet filters

### 4 pack

Use with InnoSpire Elegance Use with InnoSpire Essence Replace filter every 3 months

#### HH1312

# Maintain your InnoSpire Essence or Elegance

The air inlet filters help to keep your compressor in good working order by preventing large dust particles from entering inside. It is important to change this filter when it becomes discolored or dirty.

For maintaining your InnoSpire Essence & Elegance

Maintain your InnoSpire Essence or Elegance



## Specifications

#### **Product details**

- Air inlet filters: 4 pack
- For use with: InnoSpire Elegance, InnoSpire

Essence • Replace: 1 filter every 3 months, or when discolored or dirty

### HH1312/00

# Highlights

#### Air inlet filters - 4 pack



The pack of 4 Air inlet filters help filter large dust particles from the air that is drawn into your InnoSpire compressor, helping to keep it in good working order



Issue date 2018-07-25 © 2018 Koninklijke Philips N.V. All Rights reserved. Version: 6.1.2 Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners. 12 NC: 0000 011 02088 UPC: 3 83730 02088 4 www.philips.com