



Philips Threshold  
IMT

Inspiratory Muscle Trainer  
Increases endurance  
Increases respiratory strength

HH1332

## Pressure that's good for you

With a flow-dependent, one-way valve, threshold devices maintain a constant pressure, regardless of the patients airflow.

### Convenient design

- Reliable and easy to clean
- Works in any position

### Mouthpiece/mask option

- Mouthpiece/mask option

### Intended use

- Intended use of IMT

**PHILIPS**

IMT

Inspiratory Muscle Trainer Increases endurance, Increases respiratory strength

# Specifications

## Product details

- Warranty: Replace after 1 year
- Intended Use: As directed by your physician
- Cleaning: Wash in warm soapy water after each

- use. Flush all parts with clean water. Shake out excess water and air dry. DO NOT BOIL OR HEAT.
- Includes: training diary

HH1332/00

# Highlights

## Mouthpiece/mask option

Mouthpiece/mask option for flexibility

## Reliable and easy to clean

Threshold is made of durable, high-impact acrylic.

## Works in any position

Effective therapy in any position

## Intended use

Threshold IMT inspiratory muscle trainer can be used for muscle strength and endurance.



Issue date 2018-09-19

Version: 1.0.1

12 NC: 8670 001 46074  
EAN: 03 83730 73001 1

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)