



Philips Nuance
Mask

Standard

Pack
3 pillows (S, M, L)

A new choice for better sleep

The nasal pillow technology of Nuance delivers exceptional comfort and seal, helping you to support better sleep.



HH1100

Easy to use

- Light-weight, flexible tubing for less pull on the mask.

Appealing, minimalistic design

- Simple design is easier to use and fit.

Gel pillows

- The pillows fit difference nostril sizes.

Greater stability

- Stays in place throughout the night.

PHILIPS

Mask
Standard Pack, 3 pillows (S, M, L)

Specifications

Design

- Easy to use: Easy assembly and disassembly
- Comfort: Fabric frame

Intent of use

- To be used on patients >30kg: with CPAP prescription

Service

- Warranty: 90-days

Specifications

- Operating pressure: 4 -20 cmH2O

- Sound pressure level: 19dBA
- N° of parts: Four
- Cushion material: Silicone rubber
- Headgear material: Nylon UBL/Urethane foam/Lycra
- Frame material: Nylon, Lycra, Urethane foam and Trico mesh
- Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry
- Cleaning: headgear, tubing: Hand wash with mild detergent
- Swivel/Tubing/Elbow material: Polypropylene
- Pillows material: Urethane gel

HH1100/00

Highlights

Lightweight and simple



Easier to assemble, disassemble and use than the leading pillows mask.

Provides comfort and seal



Nuance gel pillows cushion is designed differently to provide a more comfortable experience. The gel base of the pillows should fit comfortably under your nose.

Innovative non-slip headgear



Ideal for patients that want to sleep freely.

Patient experience



During our patient preference study the Nuance platform was preferred over the leading nasal pillows mask for seal, comfort, ease of use and for less nasal irritation.



Issue date 2018-04-13

Version: 7.0.1

12 NC: 8670 001 32415

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* The purchase of this medical device does not give the right to obligatory Social Security reimbursement