

Philips
12V DC cord
DreamStation

incl. battery adapter cable
Car cable
Compatible with System One

HH1422



Ideal for frequent travelers

Part of the Dream Family

Use the cord to power your DreamStation from alternative power sources like 12 Volt DC power port in automotive vehicles or from a compatible battery. Ideal for frequent travelers, campers and long-distance drivers

The freedom to use your device on the move

- Benefit from your therapy anywhere

PHILIPS

12V DC cord DreamStation
incl. battery adapter cable Car cable, Compatible with System One

HH1422/00

Specifications

Power

- Power supply: Compatible with DreamStation and System One devices

Package contents

- 2 replacement fuses

- Battery adapter cable
- Shielded DC cord

Specifications

- Battery adapter cable: 5 m
- Length of DC cord: 2 m

Highlights

Easy to use



These cables are ideal for travelling as it gives you the option of either using a battery pack or any 12V cigarette lighter socket to power it. The alligator clips cable is for use with a deep cycle marine battery. The DC cord will power the DreamStation along with heated humidifier and heated tube if connected.



Issue date 2018-02-05

Version: 4.0.2

12 NC: 0000 011 20747
UPC: 6 06959 40756 5

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com