

Philips DreamStation Go  
Disposable 30-day fine  
filter

Filters dust, pollen, & more

Replace every 30-days  
95% efficiency at < .7 micron  
6-pack

HH1454



## 30-day fine filter

Philips 30-day disposable replacement fine filters help clean the air you breathe by capturing dust, pollen, mold spores, pet dander and more with a 95% efficiency at .5-.7 micron.

### 30-day disposable fine filter

- 30-day fine filter

# PHILIPS

Disposable 30-day fine filter  
Filters dust, pollen, & more Replace every 30-days, 95% efficiency at < .7 micron, 6-pack

HH1454/00

## Specifications

### Compatibility

- PAP System: DreamStation Go CPAP, AutoPAP

### Service

- Warranty: 15-days

### Technical information

- Replacement: Every 30 days

- Filter material: Blended Synthetic Fiber
- Filter efficiency: 95% Efficient @ .5-.7 micron

### Packaging details

- Content: Includes 2 filters

### Part Number

- Disposable 30-day fine filter: 1133744

## Highlights

### 30-day fine filter



For enhanced filtration try the 30-day disposable fine filter in your DreamStation Go travel PAP machine.



Issue date 2018-08-30

Version: 7.0.2

12 NC: 0000 011 33744  
UPC: 6 06959 41544 7

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)