




Philips DreamWear  
Full Face Magnet Clips

Full Face Replacement Part  
with 1 set of magnetic clips



HH1132

## Like wearing nothing at all.\*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

### Magnetic clips

- The headgear clips and mask cushion contain magnets

### Easy to use

- Magnet clips for headgear attachment

**PHILIPS**

Full Face Magnet Clips  
Full Face Replacement Part with 1 set of magnetic clips

HH1132/00

# Specifications

## Warning: Magnetic interference

- Beware with pacemakers: Read IFU's before use.

## Product details

- Warranty: 90 days

# Highlights

## Easy to use



The magnet clips make attaching the headgear to the frame easy and fast!

## Magnetic clips

Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants. Read in-pack IFU's for more information.



Issue date 2018-04-03

Version: 2.2.1

12 NC: 0000 011 33454  
UPC: 6 06959 41466 2

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\*\* Data analysis of Aug 2017 Patient Preference trial where n = 85 and prescribed masks include Resmed Airfit F10 and F20, Respironics Amara View and F&P Simplus; based on 10 days of use.