



Philips InnoSpire Go  
Mouthpiece assembly

Vibrating mesh technology

Replacement part  
For use with InnoSpire Go



HH1345

## Clinically proven vibrating mesh technology

The replacement mouthpiece assembly is designed to maintain the performance of InnoSpire Go and to efficiently deliver medication for children and adults.

**Easy to clean and maintain**

- Only one part to clean

**Clinically proven vibrating mesh technology**

- Vibrating mesh technology, used in hospitals worldwide

**Provides 1 year of use**

- Provides 1 year of use

**PHILIPS**

# Specifications

## Product details

- Mouthpiece assembly: for use with InnoSpire Go
- Maximum fill level: 8 ml
- Mesh material: Nickel palladium
- Outer body material: Grilamid™ (plastic)

- Packaging size: 65 mm x 65 mm x 70 mm
- Packaged weight: 45 g

## Service

- Mouthpiece assembly warranty: 1 year, recommend yearly replacement

# Highlights

## Easy to clean and maintain



The easily detachable mouthpiece makes the InnoSpire Go simple to use, clean and maintain.

## Clinically proven technology



It consistently delivers aerosol in the optimal particle size, to target the medication to your lungs.

## Therapy use



Giving you consistent and reliable therapy.



Issue date 2018-07-26

Version: 6.0.2

12 NC: 0000 011 28501  
EAN: 87 12581 75436 5

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)