



Philips DreamWear
Frame

Replacement Part
Large frame



HH1118

Like wearing nothing at all.*

With the freedom and comfort of DreamWear, you can sleep how you want.

Innovative design

- Air tube connection at the top of the head
- Soft, flexible frame

Easy to use

- Features just for you

PHILIPS

Specifications

Innovative design

- Comfort and performance: Soft fabric straps

Product details

- Materials: Polycarbonate, Silicone
- Warranty: 90 days

Cleaning: Frame


- Hand wash with mild detergent: Rinse thoroughly, air dry¹

Replacement

- Inspect daily for wear: Replace at first sign of wear

Highlights

Easy to use frame



Easy to use, choose one of the three frame size that fits you best. All the frame sizes fit all the cushion sizes.

Innovative frame design



The soft, flexible frame is comfortable to wear at night. Users have stated it feels like wearing nothing at all.*

Freedom of movement



By moving the tube location to the top of the head, DreamWear offers more freedom of movement during the night. Users stated they didn't need to choose a sleep position when wearing DreamWear.*



Issue date 2018-04-03

Version: 10.0.1

12 NC: 0000 011 16747
UPC: 6 06959 40206 5

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

** Data analysis of Aug 2017 Patient Preference trial where n = 85 and prescribed masks include Resmed Airfit F10 and F20, Respirationics Amara View and F&P Simplus; based on 10 days of use.

* *2017 Philips User Preference Questionnaire

* ¹Please refer to Instructions For Use for more detailed cleaning instructions.