

Philips System One
Ultrafine filter

6-Pack

HH1421

Additional filtration capability

Insight Driven design

Two-step air filtration system to provide optional filtration of very fine particles for your System One 60 series device.

Additional filtration of very fine particles

- Easy to insert and remove

PHILIPS

Ultrafine filter
6-Pack

Specifications

Service

- Warranty: 90 days

Design

- Ultra-fine filter: for System One 60 series

HH1421/00

Highlights

Replaceable

For optimal performance of your System One device. To be used in conjunction with the reusable pollen/foam filter to block fine particles from getting into the device. To be replaced after 30 nights of use or sooner if it shows signs of damage or appears soiled.



Issue date 2017-08-18

Version: 3.0.2

12 NC: 0000 010 29331
UPC: 6 06959 00344 6

© 2017 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com