

Replacement Part with 1 set of Fabric Wraps 2 per pack



HH1120

Like wearing nothing at all.*

With the freedom and comfort of DreamWear, you can sleep how you want.

Easy to use

· Comfortable and soft



Specifications

Product details

• Materials: Polyester, Spandex

· Warranty: 90 days

Cleaning: Fabric wraps

• Hand-wash with mild detergent: Rinse thoroughly,

hang to dry

Replacement

· Inspect daily for wear: Replace at first sign of wear

HH1120/00

Highlights

Soft fabric wraps



Provides a soft touch on your cheek



Issue date 2018-01-04

© 2018 Koninklijke Philips N.V. All Rights reserved.

Version: 7.0.1

Specifications are subject to change without notice.

12 NC: 0000 011 16754 UPC: 6 06959 02232 4

Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com

* *2017 Philips User Preference Questionnaire

* * Data analysis of Aug 2017 Patient Preference trial where n = 85
and prescribed masks include Resmed Airfit F10 and F20,
Respironics Amara View and F&P Simplus; based on 10 days of use.

* 'Please refer to the Instructions For Use for more detailed cleaning instructions.