




Philips DreamWear
Full-face cushion

Full-face replacement part
Small cushion



HH1128

Like wearing nothing at all.*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

Innovative design

- A modular design for more options

Comfort and ease of use

- Easy to use
- Under-the-nose full face cushion

PHILIPS

Full-face cushion
Full-face replacement part Small cushion

Specifications

Product details

- Materials: Polycarbonate, Silicone
- Warranty: 90 days

Cleaning: Cushion

- Hand wash with mild detergent: Rinse thoroughly,

air-dry¹

Replacement

- Inspect daily for wear: Replace at first sign of wear

HH1128/01

Highlights

Easy to use



Only a few parts for easy assembly, disassembly and cleaning.

Freedom of choice



The modular design allows users to choose the cushion that works best for them, all on one frame.

No more nose bridge red marks



The under-the-nose, innovative design of the cushion means no more red marks, discomfort and irritation on the nose bridge.**



Issue date 2018-06-22

Version: 1.2.1

12 NC: 0000 011 33430
UPC: 6 06959 41457 0

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

** Data analysis of Aug 2017 Patient Preference trial where n = 85 and prescribed masks include Resmed Airfit F10 and F20, Resprionics Amara View and F&P Simplus; based on 10 days of use.
* ¹Please refer to the Instructions For Use for more detailed cleaning instructions.