

Philips
Wisp nasal CPAP mask
with headgear

Pack
3 cushions (S/M, L, XL)
w/headgear

HH1024



Performance has a new style.

Wisp's clever auto seal groove technology accommodates your every move, for a superior seal and a great night's sleep.

Magnetic clips

- The headgear clips and mask cushion contain magnets

Open field of vision

- Read, watch TV, or wear glasses

Compact design

- Wisp comfortably fits more than 98% of users*

Performance and comfort

- Provides for quick and easy fitting and adjustment.

PHILIPS

Wisp nasal CPAP mask with headgear
Pack 3 cushions (S/M, L, XL), w/headgear

Specifications

Warning: Magnetic interference

- Beware with pacemakers: Read IFU's before use.

Compact

- Easy fitting and adjustment: Quiet exhalation

Specifications

- Sound pressure level: 19 dBA
- N° of parts: Six
- Cushion material: Silicone/polycarbonate
- Headgear material: UBL/Urethane foam/nylon/Lycra/nylon thermoplastic
- Frame material: Nylon/Spandex/Urethane foam
- Swivel material: Silicone / Polycarbonate
- Tubing material: Thermoplastic polyester elastomer
- Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry

- Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry
- Elbow/Swivel material: Polycarbonate
- Headgear clips material: Acetal copolymer
- Operating pressures: 4 -20 cmH₂O

Service

- Warranty: 90-days

Intent of use

- To be used on patients >30kg: with CPAP prescription

Design

- Comfort: Elastic crown strap
- Compact: Easy fitting and adjustment
- Patented design: Tip-of the nose cushion design
- Quiet exhalation: 360 degree rotating elbow

HH1024/00

Highlights

Greater independence



Wisp's sleek design eliminates the need for a forehead pad, offering you a better field of vision.

Compact design, natural fit



The 'tip-of-the-nose' cushion is designed to sit below your nose bridge, making minimal contact with your face while creating a superior seal. * Confirmation via the anthropometric database

Minimal design, superior seal



The design and headgear clips make it easier to use and clean.

Magnetic clips

Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants. Read in-pack IFU's for more information.



Issue date 2018-04-03

Version: 6.0.1

12 NC: 0000 010 94061
UPC: 6 06959 02127 3

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* The purchase of this medical device does not give the right to obligatory Social Security reimbursement
* Confirmation via the anthropometric database