

Philips Pico  
Nasal mask

**Pack**  
w/ Large cushion size

**HH1032**



## Value-added performance

A truly appealing traditional nasal mask that makes minimal contact with the face and offers simplicity, comfort and seal.

### **Compact design**

- Pico comfortably fits more than 98% of users\*
- The lightest and smallest traditional nasal mask \*

### **Performance and comfort**

- Provides for quick and easy fitting and adjustment

**PHILIPS**

Nasal mask  
Pack w/ Large cushion size

# Specifications

## Light and small

- Easy fitting and adjustment: Quiet exhalation

## Service

- Warranty: 90-days

## Specifications

- Operating pressure: 4 -20 cm H<sub>2</sub>O
- Sound pressure level: 19 dBA
- N° of parts: Five
- Cushion material: Silicone/polycarbonate
- Headgear material: UBL/Urethane foam/nylon/

- Lycra
- Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry
  - Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry
  - Frame/Elbow/Swivel material: Polycarbonate
  - Talon clips material: Thermoplastic polyester elastomer

## Intent of use

- To be used on patients >30kg: with CPAP prescription

HH1032/02

# Highlights

## Compact design for a natural f



The 'tip-of-the-nose' cushion is designed to sit below your nose bridge, making minimal contact with your face while creating a superior seal. \* Data on file, Internal Marketing Survey 2013

## Minimal design, superior seal



The design and headgear clips make it easier to use and clean.

## The lightest, smallest choice



Pico makes minimal contact with the face, providing the freedom to watch television or read while wearing it. \* Data on file, Internal Marketing Survey 2013



Issue date 2017-03-06

Version: 1.0.1

12 NC: 0000 011 04922  
UPC: 6 06959 00894 6

© 2017 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)