

Discover **DreamWisp**



- 1 Top-of-the-head tube connection**
Provides freedom of movement throughout the night
- 2 In-frame airflow**
Unique design allows airflow through the frame to help you sleep comfortably in any position
- 3 Four-point headgear**
Provides exceptional fit and stability
- 4 Minimal contact nasal cushion design**
For exceptional fit and comfort
- 5 Magnetic clips**
For easy assembly and release



© 2019 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.


Caution: U.S. federal law restricts these devices to sale by or on the order of a physician.

www.philips.com/respironics

RRDPGH MK 3/19/19 MCI 4108734 PN 1141861
1010 Murry Ridge Lane, Murrysville, PA 15668
800 345 6443 · 724 387 4000

PHILIPS
RESPIRONICS

DreamWisp



Comfort

with every turn

DreamWisp

Comfort with every turn

DreamWisp is our newest nasal mask specially designed with a secure, minimal contact nasal seal. This helps the mask stay securely on your face as you move during the night.

DreamWisp's design directs airflow through the frame so you can sleep comfortably, and the tubing connection located at the top of the head allows you to sleep in any position you desire — on your stomach, side or back.



See what **people**
are saying



Consumers are more satisfied with DreamWisp maintaining a seal and staying in place while sleeping.*

97% of clinicians would recommend to another clinician and a patient.**

Learn more

Contact your provider today to learn more about DreamWisp.

*Users study; Combined population of PRI Wisp and ResMed AirFit N20 users; Patient preference study, n=52, March/April 2018, at both 10 and 30 days

** April 2018 Clinician ease of use trial, n = 30; compared masks included ResMed Airfit N20 and PRI DreamWisp