



PHILIPS

Monitor/
defibrillator

Efficia DFM100

At the core of your emergency response

An emergency event is chaotic and stressful. You need to respond quickly and decisively, and your equipment needs to be ready when you are. Advanced resuscitation must be easy to learn, easy to remember, and easy to use. Philips Efficia DFM100 monitor/defibrillator tells you its readiness at a glance. Weighing only 5.66 kg (12.5 lbs), the monitor/defibrillator is lightweight enough to get to your patient quickly, and its intuitive user interface guides you through the emergency event. With the Efficia DFM100, we can help you provide dependable clinical treatment quickly and effectively.

Key advantages

- Dependable and easy to use, to be ready when you need it
- Scalable feature set to provide the level of care you need
- Enhanced cost of ownership helps you meet your budget goals with outstanding equipment

Ordering information

Part number 866199



Use model options

You must order one of these options

Options	Description
X00	<p>Hospital device Basic device, with</p> <ul style="list-style-type: none">• Manual defibrillation and synchronized cardioversion modes• 3-lead ECG monitoring• Arrhythmia detection and vitals trending <p>Includes</p> <ul style="list-style-type: none">• AC power cord• Lithium ion battery• Pads cable with test plug• Service manual CD• Therapy cable stabilizing collar• User documentation CD (with electronic Instructions for Use)• Printed quick reference card
X01	<p>Prehospital device Includes basic device (as described in option X00), and adds</p> <ul style="list-style-type: none">• Bottom cradle, connector side pouch, printer side pouch, rear pouch• Prehospital accessory storage system• Therapy cable stabilizing collar – INSTALLED

Note: C25 Docking Station Bundle must be selected for X01.

X02

Mobile device

Includes basic device (as described in option X00), and adds

- Bottom cradle, connector side pouch, printer side pouch, rear pouch
- Prehospital accessory storage system
- Therapy cable stabilizing collar—INSTALLED

Note: C25 Docking Station Bundle must be selected for X02.

Monitoring options

These options cannot be combined

Options	Description
A01	SpO₂ monitoring Pulse oximetry Includes adult reusable fingertip transducer
A02	SpO₂ and NBP monitoring Pulse oximetry and non-invasive blood pressure Includes adult re-usable fingertip transducer, adult blood pressure cuff and NBP air hose
A03	SpO₂, NBP and EtCO₂ monitoring Pulse oximetry, noninvasive blood pressure and capnography port Includes adult reusable fingertip transducer, adult blood pressure cuff and NBP air hose Note: Must be ordered with one of C07, C08, or C09*
A04	SpO₂ and EtCO₂ monitoring Pulse oximetry and capnography port Includes adult reusable fingertip transducer Must be ordered with one of C07, C08, or C09*
A05	EtCO₂ monitoring Capnography port Must be ordered with one of C07, C08, or C09*

**Note: Capnography consumables must be ordered separately.*

Therapy options

These options can be combined

Options	Description
B01	External pacing Adds noninvasive transcutaneous pacing
B02	AED mode Adds AED mode for adults and infant/child

Accessory options

These options can be combined

Options	Description
C01	External paddles Adult and infant/child external paddles with patient contact indicator (PCI), water resistant <i>Pediatric paddles are housed under the adult pads</i>
C02	USB data drive USB data drive 2.0 compatible
C03	Accessory storage system Includes monitor side storage pouch and rear storage pouch <i>Not available with option X01 or X02</i>
C04	Bedrail mount Includes one bedrail mounting bracket <i>Not available with option X01</i>
C20	Cable management straps Includes two cable straps
C22	Tray cover pouch Top pouch for prehospital carry bag
C23	Carry bag shoulder strap Shoulder strap for prehospital carry bag
C24	Tray cover pouch and shoulder strap Top pouch and shoulder strap for prehospital carry bag
C25	Docking station bundle Docking station, power cord and bedrail hook

ECG accessory options

You must order one of these options

Options	Description
C14	Trunk and 3-lead set snap ECG <i>Default option, if no other option is specified</i>
C15	Trunk and 3-lead set grabber ECG
C16	Trunk and 5-lead set snap ECG
C17	Trunk and 5-lead set grabber ECG

Capnography options

These options cannot be combined

Options	Description
C07	EtCO₂ no external module Capnography enabled; no external EtCO ₂ hardware choice <i>Can only be ordered with option A03, A04, or A05</i> <i>Not available in combination with option C08 or C09*</i>
C08	EtCO₂ external Capnostat 5 mainstream sensor Includes M2501A mainstream module <i>Can only be ordered with option A03, A04, or A05</i> <i>Not available with option C07 or C09*</i>
C09	EtCO₂ external LoFlo sidestream sensor Includes M2741A LoFlo sidestream sensor with mounting bracket <i>Can only be ordered with option A03, A04, or A05</i> <i>Not available in combination with option C07 or C08*</i>

*Note: Capnography consumables must be ordered separately.

User documentation options

These options can be combined

Options	Description
LP1	Instructions for use Includes a printed copy of the instructions for use
LP2	User training DVD DVD with self-paced instructions on how to operate the device (DVD is English only)

Service documentation options

You can only order a maximum of one of these options

Options	Description
SM1	Service manual, English Printed copy of the English service manual
SM2	Service manual, Simplified Chinese Printed copy of the Simplified Chinese service manual

Accessories

ECG trunks, lead sets and one-piece cables (3-lead)

Part number	Description
989803160641	Efficia 3/5 ECG trunk cable, AAMI/IEC
989803160651	Efficia 3-lead grabber, AAMI
989803160661	Efficia 3-lead grabber, IEC
989803170171	3-lead ECG OR trunk cable, AAMI/IEC
M1669A	3-lead ECG trunk cable, AAMI/IEC
M1671A	3-lead ICU ECG set, grabber, AAMI
M1672A	3-lead ICU ECG set grabber, IEC
M1673A	3-lead ICU ECG set, snap, AAMI
M1674A	3-lead ICU ECG set, snap, IEC
M1678A	3-lead OR ECG set, grabber, IEC

ECG trunks, lead sets and one-piece cables (5-lead)

Part number	Description
989803160691	Efficia 5-lead grabber limb, AAMI
989803160701	Efficia 5-lead grabber limb, IEC
989803170181	5-lead ECG OR trunk cable, AAMI/IEC
989803176161	5-lead snap limb, shielded, AAMI
989803176181	5-lead snap limb, shielded, IEC
M1644A	5-lead ICU ECG set, snap, limb, AAMI
M1645A	5-lead ICU ECG set, snap, limb, IEC
M1668A	5-lead ECG trunk cable, AAMI/IEC
M1968A	5-lead ICU ECG set, grabber, limb, AAMI
M1971A	5-lead ICU ECG set, grabber, limb, IEC
M1974A	5-lead OR ECG set, grabber, limb, IEC

ECG monitoring electrodes

Part number	Description
M2202A	Radiotranslucent foam monitoring electrodes (5 per pack)

Noninvasive blood pressure cuffs and Air Hoses

Part number	Description
989803177471	Reusable blood pressure interconnect tubing, 3 m (10 ft)
989803177511	Reusable blood pressure cuff - pediatric (GCF 1202)
989803177521	Reusable blood pressure cuff - adult (GCF 1203)
989803177531	Reusable blood pressure cuff - adult large (GCF 1204)
989803209761	Adult NBP Air Hose 1.5 m (new) <i>Not available for DFM100 in China</i>
989803209771	Adult NBP Air Hose 3.0 m (new) <i>Not available for DFM100 in China</i>
989803217001*	NBP inlet 5 mm (5 ea) with Allen Keys <i>Not available for DFM100 in China</i>
989803217011*	NBP inlet 5 mm (50 ea) with Allen Keys <i>Not available for DFM100 in China</i>

SpO₂ sensors and cables

Part number	Description
M1191B/ 989803144371	Reusable adult SpO ₂ finger sensor, 2 m (6.5 ft)
M1191BL 989803144381	Reusable adult SpO ₂ finger sensor, 3 m (10 ft)
M1192A/ 989803205871	Reusable pediatric/small adult finger SpO ₂ sensor, 1.5 m (5 ft)
M1196A/ 989803205901	Reusable adult SpO ₂ clip sensor, 3 m (10 ft)
M1941A/ 989803105681	SpO ₂ extension cable, 2 m (6.5 ft)

EtCO₂ mainstream sensor

Part number	Description
M2501A/ 989803142651	Capnostat 5 mainstream CO ₂ sensor

EtCO₂ mainstream filter lines

Part number	Description
M2513A/ 989803206591	Capnostat 5 reusable adult/pediatric airway adapter, ET tube size > 4.0 mm
M2533A/ 989803142661	Capnostat 5 single-patient use adult airway adapter, ET tube size > 4.0 mm

* Only used for updating to the new NBP connector and new Air Hose.

EtCO₂ sidestream sensor

Part number	Description
M2741A	LoFlo sidestream CO ₂ sensor

EtCO₂ sidestream filter lines

Part number	Description
M2744A/ 989803206611	LoFlo sidestream CO ₂ nasal cannula - adult
M2745A/ 989803206621	LoFlo sidestream CO ₂ nasal cannula - pediatric
M2750A/ 989803206641	LoFlo sidestream CO ₂ /O ₂ nasal cannula - adult
M2751A/ 989803206651	LoFlo sidestream CO ₂ /O ₂ nasal cannula - pediatric
M2756A/ 989803206671	LoFlo sidestream CO ₂ oral-nasal cannula - adult
M2757A/ 989803206681	LoFlo sidestream CO ₂ oral-nasal cannula - pediatric
M2760A/ 989803206691	LoFlo sidestream CO ₂ /O ₂ oral-nasal cannula - adult
M2761A/ 989803206701	LoFlo sidestream CO ₂ /O ₂ oral-nasal cannula - pediatric
M2768A/ 989803206711	LoFlo sidestream airway adapter set - adult/pediatric
M2776A/ 989803206751	LoFlo sidestream straight sample line

Please note not all accessories and supplies are available in all countries.
The list of accessories is subject to change without notice.

Internal paddles*

Part number	Description
M4741A	7.5 cm switched internal paddles
M4742A	6.0 cm switched internal paddles
M4743A	4.5 cm switched internal paddles

External paddles

Part number	Description
M3543A	External paddles with PCI, water resistant
M4759A*	Paddle electrode replacement for M3543A paddles
989803196431	Efficia external paddles with PCI, water resistant
989803197591	Paddle electrode replacement for 989803196431 paddles

Multi-function defibrillator pads: plug-style

Part number	Description
M3713A	Adult/child plus pads (>10 kg) (10 sets)
M3716A	Adult/child radiolucent multifunction electrode pads (10 sets)
M3717A	Infant plus pads (<10 kg) (5 sets)
M3718A*	Adult/child radiotransparent multifunction electrode pads (10 sets)
M3719A*	Infant radiotransparent multifunction electrode pads (<10 kg) (5 sets)
989803158211	HS FR/FR2 defibrillator pads (1 set)
989803158221	HS FR/FR2 defibrillator pads (5 sets)
989803139261	SMART Pads II (1 set)
989803149981	SMART Pads III (1 set)
989803149991	SMART Pads III (5 sets)

Hands-free pads therapy cables

Part number	Description
M3508A*	Hands-free cable, plug connector
989803197111	Pads adapter cable

*Not available in all markets

Cable management and accessory storage

Part number	Description
989803190331	DFM100 hospital accessory storage system; includes monitor side pouch and rear pouch
989803190341	DFM100 prehospital accessory storage system; includes therapy cable stabilizing collar, bottom cradle, monitor side pouch, printer side pouch and rear pouch
989803190351	DFM100 cable management straps; includes two straps
989803197581	DFM100 tray cover pouch
989803197571	DFM100 carry case shoulder strap
989803197671	DFM100 tray cover pouch and carry case shoulder strap
989803194281	DFM100 therapy cable collar

Data management, paper, other

Part number	Description
989803171261	USB data drive (China market only)
989803202611	USB data drive
989803190381	50 mm chemical thermal paper, no gridlines (12 rolls)

Test loads and shorted plugs

Part number	Description
M3725A	Test load, 50 ohm, plug-style
989803171271	Test plug – shorted

Power

Part number	Description
989803190371	DFM100 lithium ion battery with capacity gauge
867207	Battery charger

Mounting solutions

Part number	Description
989803190361	DFM100 bedrail mounting bracket
989803199601	DFM100 docking station bundle; includes docking station, power cord and bedrail mount (excludes Australia)
989803199661	DFM100 docking station bundle; includes docking station, power cord and bedrail mount (Australia only)
989803199241	DFM100 docking station
989803199251	DFM100 power cord for docking station (excludes Australia)
989803199611	DFM100 power cord for docking station (Australia only)
PH-0074-60	Philips HeartStart DFM100 Roll Stand <i>Not sold by Philips, please purchase Directly from GCX (www.gcx.com)</i>

Please note not all accessories and supplies are available in all countries.
The list of accessories is subject to change without notice.

Upgrades

Each upgrade can only be ordered once per device.* For each upgrade, you must order one of the available options.

Part number	Upgrade	Options	Description
866260	SpO ₂ upgrade	D01 D02 Other combinations are not possible	Adds SpO ₂ if the device has no A options Adds SpO ₂ if the device has option A05 (EtCO ₂)
866261	NBP upgrade	D01 D02 Other combinations are not possible	Adds NBP if the device has option A01 (SpO ₂) Adds NBP if the device has option A04 (SpO ₂ and EtCO ₂)
866262	EtCO ₂ upgrade	D01 D02 D03 Other combinations are not possible	Adds EtCO ₂ if the device has no A options Adds EtCO ₂ if the device has option A01 (SpO ₂) Adds EtCO ₂ if the device has option A02 (SpO ₂ and NBP)
866263	AED mode upgrade	D01 D02	Adds AED mode if the device does not have option B01 (pacing mode) Adds AED mode if the device has option B01 (pacing mode)
866264	External pacing upgrade	D01 D02	Adds pacing mode if the device does not have option B02 (AED mode) Adds pacing mode if the device has option B02 (AED mode)

* The Efficia DFM100 serial number must be provided.



Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537
 Hamburg, Germany  0123

 Philips Goldway (Shenzhen) Industrial Inc.
No.2 Keji North 3rd Road, Nanshan District Shenzhen,
P.R. China 518057

© 2023 Koninklijke Philips N.V. All rights are reserved.
Philips reserves the right to make changes in specifications and/or to discontinue
any product at any time without notice or obligation and will not be liable for any
consequences resulting from the use of this publication. Trademarks are the
property of Koninklijke Philips N.V. or their respective owners.

www.philips.com

Printed in the Netherlands.
00000219-01-00 * NOV 2023
LC2272-001-017