

# PHILIPS

IPL hair removal system

Lumea Essential

For use on body only

15 minutes to treat lower legs

Lifetime >100.000 light pulses

Extra long cord



SC1992/00



## The ideal way to prevent hair regrowth for smooth skin at home

Philips Lumea IPL works wonders to prevent the reappearance of hair on body. Gentle pulses of light, applied regularly, keep skin silky-smooth every day.

### Long lasting smoothness

- Enjoy smooth skin every day
- Effortlessly effective

### Hair removal with IPL technology

- Adapted for safe use in the comfort of your own home

### Gentle and effective, even on sensitive skin

- Five adjustable light energy settings
- For use on the body
- Bigger body attachment for fast treatment
- Adjustable to your skin type
- Slide and Flash mode for a more convenient application

### Compact corded design

- No replacement parts, no hidden costs



asimpleswitch.com

# Highlights

## Enjoy smooth skin every day



Philips Lumea applies gentle pulses of light to the hair root. This results in the hair shedding naturally and hair regrowth is inhibited. Repeating this treatment regularly will keep your skin smooth everyday.

## Effortlessly effective



Our clinical studies have shown significant hair reduction in as little as four bi-weekly treatments resulting with smooth hair free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual hair regrowth.

## IPL technology



Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop our breakthrough hair removal system. For more than 10 years we carried out extensive consumer research with more than 2000 volunteers.

## 5 light energy settings



Philips Lumea has five adjustable light energy settings to ensure gentle but effective treatment. Used correctly, Philips Lumea's IPL treatment is safe and gentle to use even on sensitive skin and sensitive body areas.

## Application areas




Use with confidence to remove hair on body areas below the neckline: legs, armpits, bikini area, arms and belly.

## Bigger body attachment, 4cm<sup>2</sup>



Bigger body attachment enables fast treatment on large areas like legs.

## Slide and Flash mode




Main application mode enables sliding motion for more convenient and continuous application.

## No replacement parts needed



Philips Lumea is a complete solution that comes ready to work straight out of the box. More importantly, it does not require any replacement lamps or gels.

## Suitable hair and skin types



Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blonde or white/grey hair. Philips Lumea is also NOT suitable for dark skin.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO<sub>2</sub> emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Technical specifications

**High performance lamp:** Delivers > 100,000 flashes

## Safety and adjustable settings

**Integrated safety system:** Prevent unintentional flashing

**5 light energy settings:** Adjustable to your skin type

**Integrated UV filter:** Protects skin from UV light

## Application areas

**Body areas:** Arms, Legs, Belly, Bikini, Armpits

## Attachments

**Body attachment (4cm<sup>2</sup>):** For use below the neckline

## Application time

**Lower leg:** 8 min

**Armpit:** 1 min

**Bikini line:** 1 min

## Items included

**Instructions for use:** Quick Start Guide, User manual

**Storage:** Pouch

**Adapter:** 24V / 1500mA

## Power

**Voltage:** 100-240 V

## Service

**Warranty:** 2 years global warranty

## Application mode

**Corded / cordless use:** Corded

**Slide and Flash:** For convenient application

**Step and Flash:** For treatment on curvy areas

