

PHILIPS

IPL hair removal
system

Lumea Precision

For use on body and face

15 minutes to treat lower legs

Lifetime >140.000 light pulses

Cordless design



SC2006/11



The ideal way to prevent hair regrowth for smooth skin at home

Philips Lumea IPL is the most effective cordless IPL solution and works wonders to prevent the reappearance of hair on body. Gentle pulses of light, applied regularly, keep skin silky-smooth every day.

Long lasting smoothness

- Enjoy smooth skin every day
- Effortlessly effective

Hair removal with IPL technology

- Adapted for safe use in the comfort of your own home
- Suitable for a wide variety of hair and skin types

Gentle and effective application on face and body

- Safe and effective even on sensitive areas
- Precision attachment for safe facial treatment
- Bigger body attachment for fast treatment
- Five adjustable light energy settings
- Slide and Flash mode for a more convenient application

Maximum freedom, minimum hassle

- Cordless convenience
- No replacement parts, no hidden costs



asimpleswitch.com

Highlights

Long lasting smoothness



Philips Lumea applies gentle pulses of light to the hair root. As a consequence the hair sheds naturally and hair regrowth is inhibited. Repeating this treatment regularly will keep your skin smooth every day.

Effortlessly effective



Our clinical studies have shown significant hair reduction in as little as four bi-weekly treatments resulting in smooth hair free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual hair regrowth.

IPL technology



Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop our breakthrough hair removal system. For more than 10 years we carried out extensive consumer research with more than 2000 women.

Suitable hair and skin types

Suitable hair and skin tones

Skin Tone	White	Beige	Light brown	Medium brown	Dark brown	Brownish black and darker
Hair Color						
Black	☑	☑	☑	☑	☑	☑
Dark brown	☑	☑	☑	☑	☑	☑
Brown	☑	☑	☑	☑	☑	☑
Dark blond	☑	☑	☑	☑	☑	☑
Light blond	☑	☑	☑	☑	☑	☑
Red	☒	☒	☒	☒	☒	☒
White/grey	☒	☒	☒	☒	☒	☒

☑ This product is suitable for you to use ☒ This product is neither suitable nor safe for your skin and hair type

Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blonde or white/grey hair. Philips Lumea is also NOT suitable for very dark skin.

Five light energy settings




Philips Lumea has five adjustable light energy settings to ensure gentle but effective treatment. Used correctly, Philips Lumea's IPL treatment is safe and gentle to use even on sensitive skin and sensitive body areas.

Application areas



Use with confidence to remove hair on face (upper lip, chin and sideburns) and body parts including legs, underarm, bikini area, belly and arms.

Cordless convenience



Lumea's unique cordless design is perfect for maximum reach and ease of use.

No replacement parts needed

Philips Lumea is a complete solution that comes ready to work right out of the box. More importantly, it does not require any replacement lamps or gels.

Extra safety for the face



The precision attachment is equipped with an extra integrated light filter and enables a safe and convenient application on the upper lip, chin and sideburns. It can also be used on other sensitive areas on the body, such as the underarms, the bikini area and hard-to-reach areas.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Technical specifications

High performance lamp: Delivers > 140,000 flashes

Safety and adjustable settings

Safety ring system: Prevent unintentional flashing

5 light energy settings: Adjustable to your skin type

Integrated UV filter: Protects skin from UV light

Application areas

Body areas: Legs, Arms, Belly, Bikini, Armpits

Face areas: Chin, Sideburns, Upper lip

Application mode

Slide and Flash: For convenient application

Step and Flash: For treatment on curvy areas

Corded / cordless use: Cordless

Attachments

Body attachment (4cm²): For use below the neckline

Precision attachment (2cm²): For safe application on face

Application time

Lower leg: 8 min

Armpit: 1 min

Bikini line: 1 min

Face areas: 1 min

Items included

Storage: Pouch

Instructions for use: Quick Start Guide, User manual

Charger: Multi voltage charger

Power

Voltage: 100-240 V

Charging time: 100 min

Running time: Minimum 370 flashes

Service

Warranty: 2 years global warranty

