

PHILIPS

Frame Connector

DreamWisp

DreamWisp Frame Replacement

Size: Medium

HH1146/01



Sleep comfortably at every turn

Combining DreamWear's top-of-the-head tube design with Wisp's minimal contact nasal cushion, DreamWisp gives you the freedom to sleep in any position throughout the night.

Minimal parts

- A desired fit with minimal parts

Innovative design

- A modular design for more options

Renew your product with original Philips parts

- Easily renew your product with original Philips parts

Specifications

Product details

Materials: Silicone

Materials: Polycarbonate

Warranty: 90 days

Calibration data

Hand wash with mild detergent: Rinse thoroughly, hang or air dry

Replacement

Inspect daily for wear: Replace at first sign of wear

