

PHILIPS

Headgear

DreamWisp

DreamWisp Headgear
Replacement

Standard size

HH1147/00



Sleep comfortably at every turn

Combining DreamWear's top-of-the-head tube design with Wisp's minimal contact nasal cushion, DreamWisp gives you the freedom to sleep in any position throughout the night.

Minimal parts


- A desired fit with minimal parts
- Easy to attach, simply to adjust

Renew your product with original Philips parts

- Easily renew your product with original Philips parts

Highlights

Designed with minimal parts



DreamWisp's head gear is designed to perform with exceptional comfort. The four-point headgear is styled to provide stability with fewer parts, while the soft, lightweight frame doubles as the DreamWisp's airflow channel. That means fewer components that come in contact with you as you sleep, so that you can enjoy comfortable therapy — all night, every night.

Replacement part

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Magnetic clips



Helping you stay committed to your sleep therapy means making it as convenient as possible. DreamWisp features magnetic clips that are fast and easy to connect or detach. That makes it simple to take your mask on or off, anytime you choose. Users stated they are most satisfied with how easy it is to use the DreamWisp mask.

Specifications

Product details

Materials: Polyurethane foam

Materials: Nylon

Materials: Spandex

Warranty: 90 days

Cleaning

Hand wash with mild detergent: Rinse thoroughly, hang or air dry

Replacement

Inspect daily for wear: Replace at first sign of wear

