

# PHILIPS

Minimal contact  
nasal mask

DreamWisp

Top-of-the-head tubing

Auto seal cushion

Easy-to-attach headgear

Available in multiple sizes

HH1142/00



## Sleep comfortably at every turn

DreamWisp was designed to help you sleep comfortably with every turn. Combining DreamWear's top-of-the-head tube design with Wisp's minimal contact nasal cushion, DreamWisp gives you the freedom to sleep in any position through the night.

### Minimal parts

- A desired fit with minimal parts
- Easy to attach, simply to adjust

### Magnetic clips

- The headgear clips and mask frame contain magnets

### Open field of vision

- Read, watch TV or wear glasses

### Freedom of movement

- Allows you to sleep in any position

### Comfort and ease of use

- The cushion with comfort and confidence

# Specifications

## Product details

**Materials:** Silicone

**Materials:** Polycarbonate

**Materials:** Polyurethane foam

**Materials:** Nylon

**Materials:** Spandex

**Materials:** Polyester

**Materials:** Magnet (Acetal)

**Warranty:** 90 days

**Operating pressure:** 4 - 30 cmH<sub>2</sub>O

**Sound level:** 18.7 dbA

## Warning: Magnetic Interference

**Beware with pacemakers:** Read IFU's before use

## Intent of use

To be used on patients (<gt/>66 lbs/30kg): for whom CPAP therapy has been prescribed

## Warning: Magnetic interference

**Beware with pacemakers:** Please refer to the IFU before use

## Cleaning

**Hand wash with mild detergent:** Rinse thoroughly, hang or air dry

## Replacement

**Inspect daily for wear:** Replace at first sign of wear

