

# PHILIPS

Minimal contact  
nasal mask

DreamWisp

Top-of-the-head tubing

Auto seal cushion

Easy-to-attach headgear

Available in multiple sizes

HH1142/00



## Sleep comfortably at every turn

DreamWisp was designed to help you sleep comfortably with every turn. Combining DreamWear's top-of-the-head tube design with Wisp's minimal contact nasal cushion, DreamWisp gives you the freedom to sleep in any position through the night.

### Minimal parts

- A desired fit with minimal parts
- Easy to attach, simply to adjust

### Magnetic clips

- The headgear clips and mask frame contain magnets

### Open field of vision

- Read, watch TV or wear glasses

### Freedom of movement

- Allows you to sleep in any position

### Comfort and ease of use

- The cushion with comfort and confidence

# Highlights

## Sleep in any position



With its revolutionary top-of-the-head mask design, the DreamWisp allows greater freedom of movement by keeping tubing out of the way. You are free to change sleeping positions throughout the night – from your side to your back to your stomach. Users stated they were satisfied with how well the mask stayed in place while sleeping.

## Designed with minimal parts



DreamWisp's head gear is designed to perform with exceptional comfort. The four-point headgear is styled to provide stability with fewer parts, while the soft, lightweight frame doubles as the DreamWisp's airflow channel. That means fewer components that come in contact with you as you sleep, so that you can enjoy comfortable therapy – all night, every night.

## Magnetic clips




Helping you stay committed to your sleep therapy means making it as convenient as possible. DreamWisp features magnetic clips that are fast and easy to connect or detach. That makes it simple to take your mask on or off, anytime you choose. Users stated they are most satisfied with how easy it is to use the DreamWisp mask.

## The Wisp cushion



DreamWisp combines the performance of a nasal mask with a minimally invasive design, an exceptional seal, and the comfort of Wisp's cushion technology. Patients are more satisfied with how well the mask maintains a seal during use.

## Unobstructed view



The unique top-of-the-head tube design allows you to have a clear field of vision. With DreamWisp, you can read, surf the web or watch TV as you wind down to sleep.

## Magnetic clips



Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants. Read in-pack IFU for more information.

# Specifications

## Product details

**Materials:** Silicone

**Materials:** Polycarbonate

**Materials:** Polyurethane foam

**Materials:** Nylon

**Materials:** Spandex

**Materials:** Polyester

**Materials:** Magnet (Acetal)

**Warranty:** 90 days

**Operating pressure:** 4 - 30 cmH2O

**Sound level:** 18.7 dbA

## Warning: Magnetic Interference

**Beware with pacemakers:** Read IFU's before use

## Intent of use

To be used on patients (<gt/>66 lbs/30kg): for whom CPAP therapy has been prescribed

## Warning: Magnetic interference

**Beware with pacemakers:** Please refer to the IFU before use

## Cleaning

**Hand wash with mild detergent:** Rinse thoroughly, hang or air dry

## Replacement

**Inspect daily for wear:** Replace at first sign of wear

