

PHILIPS

Auto CPAP

DreamStation

Colorful, easy-to-use display

Display and track nightly use

Designed to reduce noise

HH1479/00



Take charge with ease and confidence


With easy-to-navigate menus, a sleek, compact design and quiet operation, DreamStation makes it simple to start, customise and continue your obstructive sleep apnoea (OSA) therapy.

Performance

- Support for the long term
- Make progress every day
- Small and lightweight
- Performance and appeal

Highlights

Performance and appeal




DreamStation offers a sleek, low-profile design. Its vibrant color screen displays daily progress and features easy-to-use menus.

Small and lightweight




DreamStation's top-mounted humidifier hose swivel and front-facing display make it easy for you to operate while lying down or sitting up in bed. And, it's small and light, making it easy to pack for travel.

Daily Progress Feedback



Each morning you are greeted with Daily Progress Feedback – a summary of how far you've come to motivate you to stick to your therapy. The intuitive, color display shows a simple trend of their nightly hours of use, followed by a summary of the last 30 days of good nights with more than four hours of use.

Bluetooth enabled



Built-in Bluetooth makes connectivity to the DreamMapper self-management tool easy. DreamMapper provides goal-setting tools, helpful videos and ongoing feedback to keep you motivated and engaged.

Specifications

Compatibility

For use with: DreamMapper

Product details

Bluetooth: Built-in

Performance data

Flex mode: C-Flex

Ramp time: 0-45 min (5-min increments)

Starting ramp pressure: SmartRamp

Operating specifications

Altitude compensation: Automatic

Operating pressure: 4 to 20 cmH2O

