

# PHILIPS

IPL - Hair removal  
device

Lumea Prestige

For use on body and face

15 minutes to treat lower legs

>250,000 lamp flashes

SC2007/00




# Highlights

## IPL technology




Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop our breakthrough hair removal system. For more than 10 years we carried out extensive consumer research with more than 2000 volunteers.

## Enjoy smooth skin every day



Philips Lumea applies gentle pulses of light to the hair root. As a consequence, the amount of hair your body grows gradually decreases. Repeating the treatment leaves your skin beautifully hair free and touchably smooth.

## Effortlessly effective



Our clinical studies have shown significant hair reduction in as little as four fortnightly treatments resulting in smooth hair free skin. To maintain these results, simply repeat the treatment when needed. Time between

treatments may vary based on your individual hair regrowth. - To enlarge this picture, please click on the image in the image gallery on the top of this page

## Suitable hair and skin types



Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blonde or white/grey hair. Philips Lumea is also NOT suitable for very dark skin.

## Application areas



Use with confidence to remove hair on face (upper lip, chin and sideburns) and body parts including legs, underarm, bikini area, belly and arms.

## Fast treatment on legs



Big treatment window for fast application on large areas like legs.

## Extra safety for the face




The precision attachment is equipped with an extra integrated light filter and enables a safe and convenient application on the upper lip, chin and sideburns. It can also be used on other sensitive areas on the body, such as the underarms, the bikini area and hard-to-reach areas.

## Gentle treatment



Gentle to your skin.

## Cordless convenience



Lumea's unique cordless design is perfect for maximum reach and ease of use.

# Specifications

## Power

**Voltage:** 100-240 V

## Service

**Warranty:** 2 years global warranty

## Technical specifications

**High performance lamp:** Delivers <gt/> 250,000 flashes

## Application areas

**Body areas:** Arms, Legs, Belly, Bikini, Armpits

**Face areas:** Chin, Sideburns, Upper lip

## Attachments

**Precision attachment (2cm<sup>2</sup>):** For safe application on face

**Body attachment (4cm<sup>2</sup>):** For use below the neckline

## Items included

**Instructions for use:** User manual

**Charger:** Multi voltage charger

**Storage:** Luxurious pouch

**Accessory:** Cleaning cloth

## Power

**Charging time:** 100 min

**Running time:** Minimum 370 flashes

## Safety and adjustable settings

**Integrated UV filter:** Protects skin from UV light  
**5 light energy settings:** Adjustable to your skin type

**Integrated safety system:** Prevent unintentional flashing

## Application time

**Lower leg:** 8 min

**Armpit:** 1 min

**Bikini line:** 1 min

**Face areas:** 1 min

## Application mode

**Slide and Flash:** For convenient application

**Corded / cordless use:** Cordless

**Stamp and Flash:** For treatment on curvy areas

