

PHILIPS

Standard Tube

DreamStation

15 mm

Replacement tubing

HH1415



Adds to your comfort

Insight Driven design


Part of the Dream Family: designed for comfort, convenience and ease of use to enhance your therapy experience

Design

- Slim, flexible, lightweight tubing

Highlights

Freedom of movement



Not all PAP tubing is alike. Slim, lightweight and flexible 1,83m long tubing. Compatible with all DreamStation devices and fits on all masks. Non-heated tube.

Specifications

Service

Warranty: 90-days

Design

Lightweight and flexible: for DreamStation PAP devices

