



Performance has a new style.

The 'tip-of-the-nose' cushion is designed to sit below your nose bridge, making minimal contact with your face while creating an optimal seal.

Easy to use

- The headgear clips and mask cushion contain magnets
- Standard Pack fit with everything you need

Open field of vision

• Read, watch TV, or wear glasses

Compact design

• Auto seal cushion technology

Performance and comfort

• Provides for quick and easy fitting and adjustment.

On-the-nose mask HH1139/00

Highlights

Standard Pack



The Wisp Standard Pack comes with 3 cushion sizes (small/medium, large, and extra large), Fabric frame, standard headgear, and magnetic clips.

Greater independence



Wisp's sleek design eliminates the need for a forehead pad, offering you a better field of vision.

Compact design, natural fit



The 'tip-of-the-nose' cushion is designed to sit below your nose bridge, making minimal contact with your face while creating an optimal seal.

Minimal design, optimal seal



The design and headgear clips make it easier to use and clean.

Magnetic clips

Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants. Read in-pack IFU's for more information.

Specifications

Package contents

Include: 3 Nasal cushions, 1 Fabric frame, Standard headgear, Magnetic clips, User Manual

Intent of use

To be used on patients <gt/>30kg: with CPAP prescription

Cleaning: Headgear, tubing

Hand wash with mild detergent: Rinse thoroughly, hang to dry¹

Cleaning: Mask

Hand wash with mild detergent: Rinse thoroughly, air dry¹

Product details

Materials: Nylon, Polycarbonate, Silicone, Spandex, Urethane Warranty: 90 days Operating pressure: 4 to 20 cmH20

Sound level: 19 dBA Number of parts: 5

© 2023 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2023-07-03 Version: 4.4.1

UPC: 6 06959 05170 6

www.philips.com



 $^{^{}st}$ Confirmation via the anthropometric database