

PHILIPS

Disposable 30-day
fine filter

DreamStation Go

Filters dust, pollen, & more

Replace every 30-days

95% efficiency at $<lt/> .7$ micron

HH1454



30-day fine filter

Philips 30-day disposable replacement fine filters help clean the air you breathe by capturing dust, pollen, mold spores, pet dander and more with a 95% efficiency at $.5-.7$ micron.

30-day disposable fine filter

- 30-day fine filter

Disposable 30-day fine filter

HH1454/00

Highlights

30-day fine filter



For enhanced filtration try the 30-day disposable fine filter in your DreamStation Go travel PAP machine.

Specifications

Service

Warranty: 15-days

Compatibility

PAP System: DreamStation Go CPAP, AutoPAP

Technical information

Replacement: Every 30 days

Filter material: Blended Synthetic Fiber

Filter efficiency: 95% Efficient @ .5-.7 micron

Packaging details

Content: Includes 2 filters

Part Number

Disposable 30-day fine filter: 1133744

