

PHILIPS

Intense Pulsed Light
hair removal system

Lumea Comfort

For use on body

13 minutes to treat lower legs

>200,000 lamp flashes

Compact corded design + pouch

SC1985/50



Long term hair reduction

For long lasting smoothness

Philips Lumea Comfort IPL hair removal system works wonders to prevent the reappearance of hair on body. Gentle pulses of light, applied regularly, keep skin silky-smooth every day.

Long lasting hair removal with IPL technology

- 40% faster for shorter treatment time*
- Effortlessly effective
- Enjoy smooth skin every day
- Suitable for different hair and skin tones
- For effective use in the comfort of your home

Gentle and effective on body

- Gentle and effective even on sensitive areas
- Integrated Skin Tone Sensor for extra safety

Maximum freedom, minimum hassle

- Five adjustable light energy settings
- Lumea Comfort comes in a luxurious pouch
- No replacement parts, no hidden costs
- The smallest and lightest device within the Lumea range
- Slide and Flash mode for a more convenient application

Highlights

Intense Pulsed Light




Philips Lumea Comfort uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for effective use in the comfort of your own home. Philips worked closely with skin care experts to develop this technology. For more than 10 years Philips carried out extensive consumer research on IPL technology with more than 2,000 women.

Enjoy smooth skin every day



Philips Lumea Comfort applies gentle pulses of light to the hair root. This results in the hair shedding naturally and hair regrowth is inhibited. Repeating this treatment regularly will keep your skin smooth everyday.

Effortlessly effective



Our clinical studies have shown up to 75% hair reduction in as little as two months or four treatments resulting in smooth, hair-free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual

hair regrowth. - To enlarge this picture, please click on the image in the image gallery on top of this page.

40% faster*



This Lumea Comfort is 40% faster and shortens your treatment time*

Suitable hair and skin tones

Haarfarbe Haarkleur Hair color Couleur de poils	Hautton / Huidskleur / Skin tone / Couleur de la peau					
Dark blonde	✓	✓	✓	✓	✗	✗
Light blonde	✓	✓	✓	✓	✗	✗
Brown	✓	✓	✓	✓	✗	✗
Black	✓	✓	✓	✓	✗	✗
Red	✗	✗	✗	✗	✗	✗
Light blonde	✗	✗	✗	✗	✗	✗
White/grey	✗	✗	✗	✗	✗	✗

Philips Lumea Comfort works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Lumea Comfort is NOT effective on red, light-blond or white/grey hair. - To enlarge this picture, please click on the image in the image gallery on the top of this page

Application areas

Use with confidence to remove hair on parts of the body including legs, underarms, bikini area and stomach.

Skin Tone Sensor



The unique Skin Tone Sensor measures the treated skin complexion at the beginning of each session and occasionally during the

session. If it detects a skin tone that is too dark for Lumea Comfort, it will automatically stop emitting pulses.


Smallest and lightest Lumea

The Philips Lumea Comfort is the smallest and lightest in the Lumea range. It is easy to hold and easy to reach all the areas you want treat.

Luxurious pouch

This Philips Lumea Comfort has a luxurious pouch for easy storage.

Slide and Flash mode



Lumea Comfort enables a sliding motion for a more convenient and continuous application.

5 light energy settings



Philips Lumea Comfort has five adjustable light energy settings to ensure gentle but effective treatment. Used correctly, Philips Lumea's IPL treatment is comfortable and gentle to use even on sensitive skin and sensitive body areas.

No replacement parts needed

Philips Lumea Comfort is a complete solution that comes ready to work straight out of the box. More importantly, it does not require any replacement lamps or gels, providing 4 year of full body treatments - Lowest cost per treatment versus other best-selling long lasting hair removal devices.

Specifications

Application areas

Body areas: Arms, Legs, Belly, Bikini, Armpits

Safety and adjustable settings

Integrated UV filter: Protects skin from UV light

5 light energy settings: Adjustable to your skin type

Integrated safety system: Prevent unintentional flashing

Skin Tone sensor: Detects your skin tone

Application mode

Slide and Flash: For easy motion application

Corded/cordless use: Corded

Stamp and Flash: For treatment on curvy areas

Power

Voltage: 100-240 V

Technical specifications

High performance lamp: Delivers <gt/> 200,000 flashes

Service

Warranty: 2 years global warranty

Items included

Instructions for use: User manual

Storage: Luxurious pouch

Adapter: 12V / 3000mA

Application time

Lower leg: 6,5 min

Armpit: 1 min

Bikini line: 1 min

