

PHILIPS

Headgear Magnetic
Clips (2 pack)

Wisp

2 pack

HH1027



Clips make the headgear easy to adjust

Click fit headgear clips with adjustable tabs. Provides for quick and easy fitting and adjustment.

Easy to use

- The headgear clips and mask cushion contain magnets

Design

- Click fit headgear clips with adjustable tabs

Highlights

Easy to adjust



A minimal part count and click fit headgear clips make it easy to use and maintain.

Magnetic clips

Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants. Read in-pack IFU's for more information.

Specifications

Warning: Magnetic interference

Beware with pacemakers: Read IFU's before use.

Service

Warranty: 90-days

Design

Easy fitting: Self-adjusting

Specifications

Headgear clips material: Acetal copolymer and neodymium magnet

