

Philips InnoSpire Essence
Air Inlet Filters (4 pack)

4 pack

Use with InnoSpire Elegance
Use with InnoSpire Essence

HH1312



Maintain your InnoSpire Essence or Elegance

The air inlet filters help to keep your compressor in good working order by preventing large dust particles from entering inside. It is important to change this filter when it becomes discolored or dirty.

For maintaining your InnoSpire Essence & Elegance

- Maintain your InnoSpire Essence or Elegance

PHILIPS

Air Inlet Filters (4 pack)
4 pack Use with InnoSpire Elegance, Use with InnoSpire Essence

Specifications

Product details

- For use with: InnoSpire Elegance, InnoSpire Essence

- Air inlet filters: 4 pack
- Replace: 1 filter every 3 months, or when discolored or dirty

HH1312/00

Highlights

Air inlet filters - 4 pack



The pack of 4 Air inlet filters help filter large dust particles from the air that is drawn into your InnoSpire compressor, helping to keep it in good working order



Issue date 2020-06-08

Version: 6.1.2

UPC: 3 83730 02088 4

© 2020 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com