

PHILIPS

Air Inlet Filters (4 pack)

InnoSpire Essence

4 pack

Use with InnoSpire Elegance

Use with InnoSpire Essence

Replace filter every 3 months

HH1312

 Product Picture

Maintain your InnoSpire Essence or Elegance

The air inlet filters help to keep your compressor in good working order by preventing large dust particles from entering inside. It is important to change this filter when it becomes discolored or dirty.

For maintaining your InnoSpire Essence & Elegance

- Maintain your InnoSpire Essence or Elegance

Highlights

Air inlet filters - 4 pack



The pack of 4 Air inlet filters help filter large dust particles from the air that is drawn into your InnoSpire compressor, helping to keep it in good working order

Specifications

Product details

For use with: InnoSpire Elegance, InnoSpire Essence

Air inlet filters: 4 pack

Replace: 1 filter every 3 months, or when discolored or dirty

