

PHILIPS

Reusable Pollen Air
Filter

DreamStation Go

Filters dust, pollen, & more

Replace every 6 months

88% efficiency @ 7 - 10 micron

Double pack

HH1453



Reusable pollen filter

Philips reusable replacement filters help clean the air you breathe by capturing dust, pollen, mold spores, pet dander and more with an 88% efficiency at 7-10 micron.

Reusable pollen filter

- Reusable air filter

Highlights

Reusable air filter



Every 6 months replace your DreamStation Go reusable filter.

Specifications

Technical information

Replacement: Every 6 months

Filter material: 100% Polyester

Filter efficiency: 88% Efficient @ 7-10 microns

Packaging details

Content: Includes 6 filters

Part Number

Reusable pollen filter: 1133743

Service

Warranty: 90-days

Compatibility

PAP System: DreamStation Go CPAP, AutoPAP

