

# PHILIPS

Disposable 30-Day  
Fine Filter (6 pack)

DreamStation Go

Filters dust, pollen, & more

Replace every 30 days

95% efficiency at  $\leq .7$  micron

6 pack

HH1454



## 30-day fine filter

Philips 30-day disposable replacement fine filters help clean the air you breathe by capturing dust, pollen, mold spores, pet dander and more with a 95% efficiency at  $\leq .7$  micron.

### 30-day disposable fine filter

- 30-day fine filter

# Highlights

## 30-day fine filter



For enhanced filtration try the 30-day disposable fine filter in your DreamStation Go travel PAP machine.

# Specifications

## Technical information

**Replacement:** Every 30 days

**Filter material:** Blended Synthetic Fiber

**Filter efficiency:** 95% Efficient @ .5-.7 micron

## Packaging details

**Content:** Includes 2 filters

## Part Number

**Disposable 30-day fine filter:** 1133744

## Service

**Warranty:** 15-days

## Compatibility

**PAP System:** DreamStation Go CPAP, AutoPAP

