

# PHILIPS

Compressor  
Nebuliser System

InnoSpire Deluxe

Robust home use device

Treatment in under 5 minutes\*

\*Using a 2.5ml Salbutamol dose

HH1350/00



## Reliable therapy with built-in convenience

InnoSpire Deluxe is an easy to use and reliable compressor nebuliser system that includes built-in storage to keep all your therapy accessories and supplies organised and handy.

### For use in hospital and home

- Designed for continuous use

### Short treatment times

- Treatment delivered in under 5 minutes\*

### Designed for convenience

- Sturdy handle for easy transportation
- Everything you need in one place

### For commonly prescribed medications

- Compatible with commonly prescribed nebulizer medications

# Highlights

## Convenience meets reliability



InnoSpire Deluxe is designed to make your nebuliser treatment convenient and easy, with added built-in storage to house all your accessories.

## For home use



Robust enough for continuous use in both a clinical and family environment, InnoSpire Elegance provides consistent and reliable nebulizer therapy.

## Robust Design




Weighing just 4.4 lbs, InnoSpire Deluxe is light enough to pick up and move from room to room at home. With your supplies stored neatly away within the device, your nebulizer equipment stays safely together in one place.

## Treatment time



Combined with our reusable SideStream nebulizer, it is designed to provide fast drug delivery and short treatment times.

## Versatility



Capable of nebulizing a wide variety of medications that are routinely used to treat various respiratory conditions.

# Specifications

## Thermal overload protector Warranty

### Service

Warranty: 3 years

### Design

Convenient storage compartment: 2.0 kg (4.4 lbs)

### Operating parameters

Maximum flow rate: 9.8 lpm

Maximum pressure: 360 kPA (52 psi)

Electricity supply: 230V / 50Hz, 1 Amp

Flow rate: 7 lpm @ 10 psi

Thermal overload protector: Resettable thermal fuse

### Product details

Weight: 2.2 kg

Size: 380 x 170 x 110 mm

Warranty: 3 years

Noise level: 58 ±3 dBA

Treatment time (under 5 min\*): \* 2.5ml Salbutamol dose

### SideStream Plus performance

MMAD: 3.42µm

Fine particle fraction: 67%

Time (2.5mL Salbutamol): 3.68 min

### SideStream Reusable

MMAD: 2.93µm

Fine particle fraction: 76.3%

Time (2.5mL Salbutamol): 4.12 min

