

PHILIPS

Medium Headband


SmartSleep Deep
Sleep Headband

90% of users would recommend

30 day money back guarantee

Pairs with SleepMapper app

HH1600/02



Highlights

Money back guarantee

The Philips Store offers a 30 day hassle-free return policy. If you are not fully satisfied a Philips SmartSleep product within 30 days, you may contact 1-866-309-3263 for a full refund. <https://philips.com/a-w/terms-of-use.html>

SmartSleep Headband benefits



SmartSleep Deep Sleep Headband is for people who are between the ages of 18 and 50**, sleep less than 7 hours per night due to lifestyle, have no issues falling asleep or staying asleep, and sleep well when they are asleep.

SmartSleep Headband false use




SmartSleep Deep Sleep Headband does not help you fall asleep and stay asleep or prolong your deep sleep. The Deep Sleep Headband will not help with existing sleep conditions such as insomnia, restless legs or sleep apnea.

30 sensors included




SmartSleep Deep Sleep Headband uses dual sensors for highly accurate tracking that is critical to the technology. An adhesive, disposable sensor is used due to the unique placement behind the ear. Each self-adhesive sensor lasts between 1 and 3 nights of use, depending on factors such as skin conditions and time between usages, when properly applied and stored. Experience the sleep sensing technology of SmartSleep Deep Sleep Headband. The Deep Sleep Headband comes with 30 sensors, replenishments sold separately.

Clinically proven



SmartSleep Deep Sleep Headband is clinically proven to help users improve their sleep quality, increase energy, boost alertness, and reduce daytime sleepiness. 70% of users reported feeling less tired during the day.*

SleepMapper mobile app



The SleepMapper mobile app tracks your sleep patterns over time, measures key metrics and shows the slow wave boost you get from SmartSleep Deep Sleep Headband each night. Sync the app to your headband each morning to see the enhancement given each night. Please consult the Google Play Store or Apple App Store to check the compatibility of your operating system before purchase.

Adjustable sliders



Approximately 70% of users choose size Medium headband, while 30% choose a Large headband. Size Medium headband fits a head circumference of 21-23.5 inches (54-60cm). Size Large headband fits a head circumference of 23.5-26 inches (60-66cm).

Specifications

Technical specifications

Voltage: 120V AC

Technical specifications

Warranty: 2 years

Package contents

Include: SmartSleep headband, 90-days of SmartSleep Sensors, Charger, Quickstart Guide, User Manual

Material

Made with: Silver

Not made with: DEHP or latex



- * When used for 2 weeks in people who do not get enough sleep due to lifestyle
- ** Note: Sensing accuracy and product performance is lower in people over age 50 due to natural reduction in slow wave amplitude. Volume does not exceed 80dB. People with hearing impairment unable to hear 80dB will not benefit from SmartSleep Deep Sleep Headband.
- * SmartSleep Deep Sleep Headband is for people who rarely get 6 hours of sleep or people who consistently sleep 5 to 7 hours or do not get enough sleep due to lifestyle.