

# PHILIPS

Headgear (L)

Wisp

Large

HH1006



## Click fit headgear for easy adjustment

Elastic crown strap with integrated tube holder. Self adjusting for a more natural fit.

### Design

- Personalized comfort fit

# Highlights

## Easy to adjust

Self-adjusting top strap with tube holder loop

# Specifications

## Service

**Warranty:** 90-days

## Design

**Comfort:** Self-adjusting

## Specifications

**Headgear material:** UBL/Urethane foam/nylon/Lycra/nylon thermoplastic

**Cleaning: headgear, tubing:** Hand wash with mild detergent, Rinse thoroughly, hang to dry

