

PHILIPS

Headgear
(Standard)

Amara View

Standard size

HH1230



Enjoy the view

Full face mask offering comfort, performance and an open field of vision

Easy to use

- Comfort and performance

Highlights

Innovative design



One modular frame for all cushion sizes

Specifications

Innovative design

Comfort and performance: Soft fabric straps

Specifications

Headgear material: Polyurethane, Nylon-Lycra

Cleaning headgear: Hand wash with mild detergent, Rinse thoroughly, hang to dry

Service

Warranty: 90-days

