

# PHILIPS

Minimal Contact Nasal  
Mask - Fit Pack (P, S/M, L)

Wisp

Fit Pack

3 cushions (P, S/M, L)

HH1037



## Performance has a new style.

Wisp's clever auto seal groove technology accommodates your every move, for a superior seal and a great night's sleep.

### **Open field of vision**

- Read, watch TV, or wear glasses

### **Compact design**

- Wisp comfortably fits more than 98% of users\*

### **Performance and comfort**

- Provides for quick and easy fitting and adjustment.


# Highlights

## Compact design, natural fit



The 'tip-of-the-nose' cushion is designed to sit below your nose bridge, making minimal contact with your face while creating a superior seal. \* Confirmation via the anthropometric database

## Greater independence



Wisp's sleek design eliminates the need for a forehead pad, offering you a better field of vision.

## Minimal design, superior seal



The design and headgear clips make it easier to use and clean.

# Specifications

## Service

Warranty: 90-days

## Design

Comfort: Elastic crown strap

## Specifications

Operating pressure: 4 -30 cmH<sub>2</sub>O

Sound pressure level: 19 dBA

N° of parts: Six

Cushion material: Silicone/polycarbonate

Headgear material: UBL/Urethane

foam/nylon/Spandex/nylon thermoplastic

Tubing material: Thermoplastic polyester elastomer

Frame material: Nylon/Spandex/Urethane foam

Swivel material: Silicone / Polycarbonate

Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry

Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry

Elbow/Swivel material::

Silicone/Polycarbonate

Operating pressures: 4 -20 cmH<sub>2</sub>O

Headgear clips material: Acetal copolymer

## Intent of use

To be used on patients <gt/>30kg: with CPAP prescription

## Compact

Easy fitting and adjustment: Quiet exhalation

