

# PHILIPS

Overnight Battery

DreamStation Go

Maintain therapy off the grid

Compact

Lightweight, 1.53 lbs

HH1450



## CPAP therapy even 'off the grid'

For many sleep apnea patients, travelling off the power grid can create problems with maintaining therapy. Philips developed the DreamStation Go overnight battery.

**Compact and lightweight**

- Compact and lightweight

## Highlights

### Compact and lightweight

If you plan on traveling off the power grid, DreamStation Go's overnight battery can help you stay committed to your therapy. DreamStation Go set at 10cm of pressure displayed an average run time of 13 hours\* when running off the Dreamstation Go overnight battery. Note: Battery run times may

vary due to higher CPAP pressure settings, excessive mask leakage, altitude, temperature, charging devices with USB port and others. \*PAP run-time test conditions: CPAP mode, pressure 10 cm H<sub>2</sub>O, 12mm tubing, 37LPM of leak, 73.4 F room at an elevation of 1650 ft.

## Specifications

### Technical information

**Weight:** 696g (1.53lb)  
**Output voltage:** 18-24.6 VDC  
**Battery technology:** Lithium Ions  
**Input Voltage Range:** 100-240 VAC, 50/60Hz, 2.0-1.0A  
**Output power:** 50W  
**Run time:** > 8 hours  
**Minimum life cycle:** 500 cycles  
**Recharge time:** < /> 5 hours

### Part Number

**Battery Pack:** 1133281

### Service

**Warranty:** 1 year

### Compatibility

**PAP System:** DreamStation Go CPAP, AutoPAP

