

# PHILIPS

Tube & Elbow  
Assembly

Wisp

Standard

HH1008



## Connects to the front of the mask

Connect the Wisp tubing and elbow swivel to the front of the mask. You may also slide the mask tubing through the headgear crown loop before pulling on the headgear

### Design

- Rotating elbow

# Highlights

## Rotating elbow

360° rotating elbow features a quick release design with quiet exhalation

# Specifications

## Service

Warranty: 90-days

## Design

Quick release: 360° rotating elbow

## Specifications

Tubing material: Thermoplastic polyester elastomer

Elbow/Swivel material:: Polycarbonate

