

PHILIPS

Over-the-nose
nasal mask

Pico

Traditional nasal mask

Tube-in-front design

Adjust with headgear clips

HH1032/01



Value-added performance

A truly appealing traditional nasal mask that makes minimal contact with the face and offers simplicity, comfort and seal.

Compact design

- Pico comfortably fits more than 98% of users*
- The lightest and smallest traditional nasal mask *

Performance and comfort

- Provides for quick and easy fitting and adjustment

Highlights

Compact design for a natural fit




The 'tip-of-the-nose' cushion is designed to sit below your nose bridge, making minimal contact with your face while creating a superior seal. * Data on file, Internal Marketing Survey 2013

Minimal design, superior seal



The design and headgear clips make it easier to use and clean.

The lightest, smallest choice



Pico makes minimal contact with the face, providing the freedom to watch television or read while wearing it. * Data on file, Internal Marketing Survey 2013

Specifications

Service

Warranty: 90-days

Specifications

Operating pressure: 4 -20 cm H2O

Sound pressure level: 19 dBA

N° of parts: Five

Cushion material: Silicone/polycarbonate

Headgear material: UBL/Urethane foam/nylon/Lycra

Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry

Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry

Frame/Elbow/Swivel material: Polycarbonate

Talon clips material: Thermoplastic polyester elastomer

Intent of use

To be used on patients <gt/>30kg: with CPAP prescription

Light and small

Easy fitting and adjustment: Quiet exhalation

