

PHILIPS

Cushion (S/M)

Wisp

Small/Medium

HH1003



Soft and comfortable for restful night's sleep

Patented 'Tip-of the nose' cushion design, featuring auto seal groove technology. Delivers a minimal contact experience with an exceptional seal.

Design

- Less nostril irritation

Cushion (S/M)

HH1003/00

Highlights

Soft and comfortable

All three cushion sizes fit on the same frame

Specifications

Service

Warranty: 90-days

Design

Size: Small/medium

Specifications

Cushion material: Silicone/Polycarbonate

© 2021 Koninklijke Philips N.V.
All Rights reserved.

Issue date 2021-04-29
Version: 3.0.3

Specifications are subject to change
without notice. Trademarks are the
property of Koninklijke Philips N.V. or
their respective owners.

EAN: 00 60695 94074 28

www.philips.com

