

PHILIPS

Nasal Mask - Fit
Pack (S/M, L, XL)

Wisp

3 cushions (S/M, L, XL)

With silicone frame

HH1023



Performance has a new style

Wisp's clever auto seal groove technology accommodates your every move, for a superior seal and a great night's sleep.

Magnetic clips

- The headgear clips and mask cushion contain magnets

Open field of vision

- Read, watch TV, or wear glasses

Compact design

- Wisp comfortably fits more than 98% of users*

Performance and comfort

- Provides for quick and easy fitting and adjustment.


Highlights

Compact design, natural fit



The 'tip-of-the-nose' cushion is designed to sit below your nose bridge, making minimal contact with your face while creating a superior seal. * Confirmation via the anthropometric database

Greater independence



Wisp's sleek design eliminates the need for a forehead pad, offering you a better field of vision.

Minimal design, superior seal



The design and headgear clips make it easier to use and clean.

Magnetic clips

Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants. Read in-pack IFU's for more information.

Specifications

Warning: Magnetic interference

Beware with pacemakers: Read IFU's before use.

Service

Warranty: 90-days

Design

Compact: Easy fitting and adjustment

Comfort: Elastic crown strap

Patented design: Tip-of the nose cushion design

Quiet exhalation: 360 degree rotating elbow

Silicone frame: Clear unobtrusive frame

Specifications

Sound pressure level: 19 dBA

N° of parts: Six

Cushion material: Silicone/polycarbonate

Headgear material: UBL/Urethane foam/nylon/Spandex/nylon thermoplastic

Tubing material: Thermoplastic polyester elastomer

Frame material: Silicone/Polycarbonate

Swivel material: Silicone / Polycarbonate

Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry

Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry

Elbow/Swivel material:: Polycarbonate

Operating pressures: 4 -20 cmH2O

Headgear clips material: Acetal copolymer

Intent of use

To be used on patients <gt/>30kg: with CPAP prescription

Compact

Easy fitting and adjustment: Quiet exhalation

