

# PHILIPS

Reusable pollen air  
filter

DreamStation Go

Filters dust, pollen, & more

Replace every 6 months

88% efficiency @ 7 - 10 micron

Double pack

HH1453



## Reusable pollen filter

Philips reusable replacement filters help clean the air you breathe by capturing dust, pollen, mold spores, pet dander and more with an 88% efficiency at 7-10 micron.

### **Reusable pollen filter**

- Reusable air filter

# Highlights

## Reusable air filter



Every 6 months replace your DreamStation Go reusable filter.

# Specifications

## Service

Warranty: 90-days

## Compatibility

PAP System: DreamStation Go CPAP, AutoPAP

## Technical information

Replacement: Every 6 months

Filter material: 100% Polyester

Filter efficiency: 88% Efficient @ 7-10 microns

## Packaging details

Content: Includes 6 filters

## Part Number

Reusable pollen filter: 1133743

