



PHILIPS

Trilogy EV300
ventilator

Ordering and
accessories guide



Contents

Ventilator unit	6
Power cords	8
O ₂ hoses	9
Instructions for use	10
Roll stand and accessories	11
Patient circuits	14
Humidification devices	20
Patient interfaces	22
Additional parts and accessories	36
Part number index	39
Alphabetic index	42



Trilogy EV300 ventilator

The Trilogy EV300 ventilator provides invasive and noninvasive, positive-pressure ventilation to pediatric through adult patients with a minimum weight of 2.5 kg. It is an electronically controlled, pneumatic ventilation system with an integrated air compressing system. It is compatible with a range of accessories to provide a variety of therapy modes.



Features

Ventilation

- Invasive and noninvasive volume and pressure modes
- Leak compensation through Auto-Trak, which automatically adjusts trigger and cycle thresholds
- Compatible with patient circuits with passive and active exhalation valves
- Suitable for adult and pediatric patients down to 2.5 kg

Monitoring

- Numeric patient data display
- Digital airway pressure manometer
- Direct, proximal flow and volume measurements when used with active flow circuit
- FiO₂ monitoring
- EtCO₂ with CO₂ accessory
- SpO₂ and pulse rate with pulse oximeter accessory

Electrical and gas supplies

- Internal and external, hot-swappable lithium-ion batteries are standard features
- Integral O₂ blender uses either high-pressure or low-flow O₂ input
- Integral blower

Mounting

- Roll stand with O₂ cylinder holder and O₂ manifold

Accessories

- Compatible with dual- and single-limb circuits
- Compatible with passive circuits, which incorporate exhalation ports, and active circuits, which use exhalation valves

Others

- Compact and lightweight with integrated carrying handle for transport
- Identical user interface, modes, and circuit configurations as its homecare namesake, to facilitate the hospital to home transition
- Remote alarm connectivity

Configuring the Trilogy EV300

Item	Description	Page
1	Ventilator unit	5
2	Power cord*	8
3	O ₂ hose and fitting*	9
4	Instructions for use*	10
5	Roll stand and accessories	11
6	Patient circuit	14
7	Humidification device	20
8	Patient interface	22
9	Additional parts and accessories	36


* Included in region-specific ventilator units. If you order the international ventilator unit, order these parts separately.

1 Ventilator unit

Order an international version unless a region-specific version is listed. The contents of the ventilator unit are listed on p. 5.

Description	Order number	Part number
Trilogy EV300 ventilator, international	989805661891	IN2200X15B
Trilogy EV300 ventilator, international*	989805667221	FX2200X15B
Trilogy EV300 ventilator, Benelux	989805663381	BL2200X15B
Trilogy EV300 ventilator, Brazil	989805661841	BR2200X18B
Trilogy EV300 ventilator, Canada	989805661851	CA2200X12B
Trilogy EV300 ventilator, China	989805661861	CN2200X17B
Trilogy EV300 ventilator, Eastern Europe	989805663381	EE2200X15B
Trilogy EV300 ventilator, EU	989805663401	EU2200X15B
Trilogy EV300 ventilator, EU (non-Bluetooth)	989805663411	EU2200X19
Trilogy EV300 ventilator, France	989805663421	FR2200X14B
Trilogy EV300 ventilator, Germany	989805661871	DE2200X13B
Trilogy EV300 ventilator, Great Britain	989805663431	GB2200X15B
Trilogy EV300 ventilator, Iberia	989805663391	ES2200X15B
Trilogy EV300 ventilator, India	989805663441	IA2200X15B
Trilogy EV300 ventilator, Italy	989805663451	IT2200X21B
Trilogy EV300 ventilator, Korea	989805661911	KR2200X15B
Trilogy EV300 ventilator, Nordics	989805663461	ND2200X15B
Trilogy EV300 ventilator, Russia	989805663001	RU2200X15B
Trilogy EV300 ventilator, Turkey	989805663011	TR2200X15B
Trilogy EV300 ventilator, USA	989805661881	DS2200X11B

* Country exclusions may apply.




Contents of ventilator unit

International device

The international Trilogy EV300 includes the following.

Description	Order number	Part number
Trilogy EV300 ventilator unit. A DISS male O ₂ fitting comes standard on the ventilator unit. Order an additional, region-specific fitting if needed (p. 9).	—	—
Filter, particulate, each	Part of 989805662301	Part of 1134430
Patient circuit, single-limb, passive, non-heated, single-use, 22-mm ID, each	Part of 989805662101	Part of 1132340
Clamp, power cord (Not shown)	989805662271	1134347




Contents of international EV300

Region-specific device

The region-specific Trilogy EV300 includes the following.

Description	Order number	Part number
Trilogy EV300 ventilator unit	—	—
Filter, particulate, each	Part of 989805662301	Part of 1134430
Patient circuit, single-limb, passive, non-heated, single-use, 22-mm ID, each	Part of 989805662101	Part of 1132340
Clamp, power cord (Not shown)	989805662271	1134347
Power cord	See table on p. 7	
O ₂ hose		
O ₂ fitting		
Instructions for use		
O ₂ input quick connect (Not shown)		
Flex tracheal adapter (Not shown)		



Contents of region-specific Trilogy EV300. Unlabeled parts vary by region. Some versions contain additional parts.









Region-specific ventilator contents

Region	Power cord Order number (Part number) (p. 8)	O ₂ hose Order number (Part number) (p. 9)	O ₂ fitting Order number (Part number) (p. 9)	Instructions for use Order number/ (Part number) (p. 10)	O ₂ input quick- connect (p. 37)	Flex tracheal adapter (p. 17)
Benelux	989805662171 (1133284)	989805613851 (1054836)	989805621691 (1071790)	989805662661 (1143299) 989805662551 (1143288) 989805662491 (1143282)	✓	—
Brazil	989805662211 (1133468)	989805616881 (1001664)	—	989805662701 (1143303)	—	—
Canada	989805662201 (1133287)	989805616201 (1062978)	989805621701 (1071791)	989805662471 (1143280) 989805662561 (1143289)	—	✓
China	989805662221 (1133469)	989805613851 (1054836)	989805621691 (1071790)	989805662781 (1143311)	—	—
Eastern Europe	989805662171 (1133284)	989805613851 (1054836)	989805621691 (1071790)	989805662441 (1143217) 989805662711 (1143304) 989805662731 (1143306)	✓	—
EU	989805662171 (1133284)	989805613851 (1054836)	989805621691 (1071790)	—	✓	—
France	989805662171 (1133284)	989805613851 (1054836)	989805621691 (1071790)	989805662551 (1143288)	✓	—
Germany	989805662171 (1133284)	989805612551 (1054833)	989805621691 (1071790)	989805662491 (1143282)	✓	—
Great Britain	989805662241 (1133471)	989805615961 (1054832)	989805621691 (1071790)	989805662471 (1143280)	✓	—
Iberia	989805662171 (1133284)	989805613851 (1054836)	989805621691 (1071790)	989805662521 (1143285) 989805662691 (1143302)	✓	—
India	989805662251 (1133472)	989805613851 (1054836)	989805621691 (1071790)	989805662471 (1143280)	—	—
Italy	989805662171 (1133284)	989805615971 (1054834)	989805621691 (1071790)	989805662601 (1143293)	✓	—
Korea	989805662171 (1133284)	989805616881 (1001664)	—	989805662621 (1143295)	—	—
Nordics	989805662171 (1133284)	989805615961 (1054832)	989805621691 (1071790)	989805662451 (1143278) 989805662541 (1143287) 989805662671 (1143300) 989805662761 (1143309)	✓	—
Russia	989805662171 (1133284)	989805616881 (1001664)	—	989805662721 (1143305)	✓	—
Turkey	989805662241 (1133471)	989805616881 (1001664)	—	989805662771 (1143310)	✓	—
USA	989805662201 (1133287)	989805616881 (1001664)	—	989805662461 (1143279)	—	✓

2 Power cords

Power cords are 2.44 m (8 ft) long unless otherwise indicated.


Region	Description	Order number	Part number
Argentina	Power cord, Trilogy EV300, Argentina	989805662381	1136845
Australia, New Zealand	Power cord, Trilogy EV300, Australia	989805662191	1133286
Brazil	Power cord, Trilogy EV300, Brazil	989805662211	1133468
China, Hong Kong, Taiwan	Power cord, Trilogy EV300, China	989805662221	1133469
Armenia, Austria, Azerbaijan, Belgium, Bolivia, Bosnia-Herzegovina, Bulgaria, Chile, Croatia, Czech Republic, Denmark, Egypt, Estonia, Finland, France, Germany, Greece, Hungary, Indonesia, Italy, Latvia, Libya, Lithuania, Macedonia, Mongolia, Morocco, Netherlands, Norway, Peru, Poland, Portugal, Romania, Russia, Serbia, Slovakia, Slovenia, South Korea, Spain, Sweden, Switzerland, Tajikistan, Tunisia, Turkey, Turkmenistan, Vietnam	Power cord, Trilogy EV300, continental Europe (Schuko)	989805662171	1133284
India	Power cord, Trilogy EV300, India	989805662251	1133472
USA, Canada, Colombia, Honduras, Mexico, Philippines, Puerto Rico, Thailand	Power cord, Trilogy EV300, North America	989805662201	1133287
	Power cord, Trilogy EV300, North America, 4-m (13-ft)	989805662231	1133470
Bahrain, Ireland, Kenya, Kuwait, Lebanon, Oman, Malaysia, Qatar, Saudi Arabia, Singapore, United Arab Emirates, United Kingdom	Power cord, Trilogy EV300, United Kingdom	989805662241	1133471



3 O₂ hoses


Hoses are 2.5 m (8 ft) long unless otherwise indicated.

Region	Description	Order number	Part number
Australia, New Zealand	Hose, O ₂ , SIS x SIS, white	989805616871	1001659
Denmark, Finland, Norway, Sweden	Hose, O ₂ , NIST x AGA, white	989805615961	1054832
Canada	Hose, O ₂ , DISS male x DISS female, 3-m (10-ft)	989805616201	1062978
Singapore	Hose, O ₂ , DISS female x DISS female, 3-m (10-ft)	989805612421	1002556
Armenia, Austria, Azerbaijan, Germany, Libya, Mongolia, Russia, Switzerland, Tajikistan, Turkmenistan	Hose, O ₂ , NIST x DIN, white	989805612551	1054833
Bahrain, Egypt, Ireland, Kenya, Kuwait, Lebanon, Oman, Qatar, Saudi Arabia, United Arab Emirates, United Kingdom	Hose, O ₂ , NIST x BS, white	989805612691	1001663
Greece, Italy, Morocco, Tunisia	Hose, O ₂ , NIST x AFNOR, white	989805612701	1054835
Italy, Portugal, Spain	Hose, O ₂ , NIST x UNIFOR, white	989805615971	1054834
Argentina, Bolivia, Brazil, Chile, Colombia, Honduras, Indonesia, Israel, Mexico, Panama, Peru, Philippines, Puerto Rico, South Korea, Thailand, Turkey, USA, Vietnam	Hose, O ₂ , DISS female x DISS female, green, 3-m (10-ft)	989805616881	1001664
Continental Europe, China, Hong Kong, India, Ireland, Malaysia, Pakistan, South Africa, Taiwan	Hose, O ₂ , NIST x NIST, white	989805613851	1054836



O₂ fittings

Description	Order number	Part number
Fitting, O ₂ , Trilogy, DISS male	989805621681	1071789
Fitting, O ₂ , Trilogy, DISS female (Canada)	989805621701	1071791
Fitting, O ₂ , Trilogy, NIST	989805621691	1071790
Fitting O ₂ , Trilogy, SIS	989805621711	1071792





4 Instructions for use

Additionally, electronic copies are available at www.philips.com/hrcmanuals.

Description	Order number	Part number
Instructions for use, Trilogy EV300, English, non-U.S.	989805662471	1143280
Instructions for use, Trilogy EV300, English, U.S.	989805662461	1143279
Instructions for use, Trilogy EV300, Arabic	989805662481	1143281
Instructions for use, Trilogy EV300, Bulgarian	989805662431	1143216
Instructions for use, Trilogy EV300, simplified Chinese	989805662781	1143311
Instructions for use, Trilogy EV300, traditional Chinese	989805662791	1143312
Instructions for use, Trilogy EV300, Croatian	989805662581	1143291
Instructions for use, Trilogy EV300, Czech	989805662441	1143217
Instructions for use, Trilogy EV300, Danish	989805662451	1143278
Instructions for use, Trilogy EV300, Dutch	989805662661	1143299
Instructions for use, Trilogy EV300, Estonian	989805662531	1143286
Instructions for use, Trilogy EV300, Finnish	989805662541	1143287
Instructions for use, Trilogy EV300, French	989805662551	1143288
Instructions for use, Trilogy EV300, French (Canada)	989805662561	1143289

Description	Order number	Part number
Instructions for use, Trilogy EV300, German	989805662491	1143282
Instructions for use, Trilogy EV300, Greek	989805662511	1143284
Instructions for use, Trilogy EV300, Hebrew	989805662571	1143290
Instructions for use, Trilogy EV300, Hungarian	989805662591	1143292
Instructions for use, Trilogy EV300, Italian	989805662601	1143293
Instructions for use, Trilogy EV300, Korean	989805662621	1143295
Instructions for use, Trilogy EV300, Latvian	989805662641	1143297
Instructions for use, Trilogy EV300, Lithuanian	989805662631	1143296
Instructions for use, Trilogy EV300, Macedonian	989805662651	1143298
Instructions for use, Trilogy EV300, Norwegian	989805662671	1143300
Instructions for use, Trilogy EV300, Polish	989805662681	1143301
Instructions for use, Trilogy EV300, Portuguese	989805662691	1143302
Instructions for use, Trilogy EV300, Brazilian Portuguese	989805662701	1143303
Instructions for use, Trilogy EV300, Romanian	989805662711	1143304
Instructions for use, Trilogy EV300, Russian	989805662721	1143305
Instructions for use, Trilogy EV300, Serbian	989805662751	1143308
Instructions for use, Trilogy EV300, Slovak	989805662731	1143306
Instructions for use, Trilogy EV300, Slovenian	989805662741	1143307
Instructions for use, Trilogy EV300, Spanish	989805662521	1143285
Instructions for use, Trilogy EV300, Swedish	989805662761	1143309
Instructions for use, Trilogy EV300, Turkish	989805662771	1143310

5 Roll stand and accessories

Order a roll stand, and configure it with any of the parts that follow.

Roll stand

Includes a basket, a slide bracket mount for a humidifier, cord and hose wraps, side mounting rails, and mounting plate. The roll stand comes partially assembled.

Description	Order number	Part number
Roll stand, small product, EV300, 111-cm (43.7-in.) high	989805660781	1141765
Roll stand, small product, EV300, 104-cm (41-in.) high	989805660791	1141766

O₂ cylinder holder

Mounts two O₂ cylinders to the roll stand. Also order the O₂ manifold/transport kit to allow seamless transitions between wall O₂ and cylinder O₂.

Description	Order number	Part number
O ₂ cylinder holder, small product roll stand, for E/M24 size cylinders (11-cm (4.3-in.) diameter)	989805656081	1131463
O ₂ cylinder holder, small product roll stand, for 14-cm (5.5-in.) diameter cylinders	989805656061	1132041

O₂ manifold/transport kit


Allows easy transition between wall O₂ and two O₂ cylinders without interruption. Kits include cylinder connecting hoses. Also order the O₂ cylinder holder to mount the cylinders to the roll stand.

Description	Order number	Part number
O ₂ manifold/transport kit, for DISS male O ₂ fitting, with green hoses	989805648241	1109602
O ₂ manifold/transport kit, for NIST O ₂ fitting, with white hoses	989805648741	1109881
O ₂ manifold/transport kit, for SIS O ₂ fitting, with white hoses	989805650181	1113371
O ₂ manifold/transport kit, for DISS female O ₂ fitting, with white hoses (Canada)	989805650191	1113392

Small products roll stand conversion kit

This kit permits the small product roll stand, PN 1126379 or 1131462, to accommodate the Trilogy EV300.

Description	Order number	Part number
Conversion kit, small products roll stand to Trilogy EV300	989805660771	1141764



Support arm

The support arm screws into the roll stand.

Description	Order number	Part number
Support arm, three-section	989805617961	332497



Tubing hangers

Description	Order number	Part number
Tubing hanger, single-limb, 22-mm ID tube	989805648891	311107
Tubing hanger, dual-limb, 22-mm ID tubes	453561507401	1003661


CleanClip mask storage system

Used to store an AF541, AF531, or PerforMax pediatric mask (i.e., one that includes a plastic shell). Attaches to either a rail or humidifier slide bracket.

Description	Order number	Part number
Bracket, CleanClip mask storage	989805627181	1076430



CleanClip mask storage system



Fisher & Paykel humidifier/water bag holder

Mounts the water bag holder and humidifier to the humidifier slide bracket of a roll stand.

Description	Order number	Part number
Bracket, water bag holder mounting, pole-mount	989805613131	1002229
Water bag holder, bent	989805613541	1002506

Mounts the water bag holder to the rail of a roll stand using the accessory post.

Description	Order number	Part number
Post, accessory, rail-mount	453561506591	1048899
Mounting kit, Fisher & Paykel humidifier with water bag holder. Includes:	989805613291	1009753
• Clamp, pole-mount	989805613311	1002230
• Bracket, water bag holder mounting, pole-mount	989805613131	1002229
• Water bag holder, bent	989805613541	1002506


Mounts the water bag holder to the stand rail using the support arm bracket.

Description	Order number	Part number
Bracket, support arm, rail-mount	989805611511	1002497
Water bag holder, bent	989805613541	1002506


Water bag holder, pole-mount

Mounts the water bag holder to the humidifier slide bracket of a roll stand using the accessory post.

Description	Order number	Part number
Bracket, water bag holder mounting, pole-mount, C-clamp style	989805605731	1052863



Water bag holder and humidifier mounted to stand with pole-mount bracket



Mounted to roll stand rail using accessory post

6 Patient circuits


Passive single-limb circuits

Passive single-limb circuits include an exhalation port. In all passive circuits, the exhalation port is a swivel port.

Heated circuits for use with WILAmed AIRcon humidifier


WILAmed heated circuits are not available in North America.

Description	Quantity	Order number	Part number
Patient circuit, single-limb, passive, heated, single-patient use, 22-mm ID (WILAmed)	10	989805662111	1132341
Patient circuit, single-limb, passive, heated, single-patient use, 15-mm ID (WILAmed)	10	989805662121	1132342



Non-heated circuits

Description	Quantity	Order number	Part number
Patient circuit, single-limb, passive, non-heated, single-patient use, 15-mm ID	10	989805662091	1132339
Patient circuit, single-limb, passive, non-heated, single-patient use, 22-mm ID	10	989805662101	1132340



Exhalation ports

For use with all passive circuits. Not available in North America.

Description	Quantity	Order number	Part number
Exhalation port, swivel, disposable (WILAsilent)	10	989805643161	1139909





Active single-limb circuits

Active single-limb circuits include an exhalation valve.

Heated circuits for use with WILAmed AIRcon humidifier



WILAmed heated circuits are not available in North America.

Description	Quantity	Order number	Part number
Patient circuit, single-limb, with active exhalation valve, heated, single-use, 22-mm ID (WILAmed)	10	989805662161	1132346
Patient circuit, single-limb, with active exhalation valve, heated, single-use, 15-mm ID (WILAmed)	10	989805662151	1132345



Non-heated circuits

Description	Quantity	Order number	Part number
Patient circuit, single-limb, with active exhalation valve, non-heated, single-use, 22-mm ID	10	989805662141	1132344
Patient circuit, single-limb, with active exhalation valve, non-heated, single-use, 15-mm ID	10	989805662131	1132343





Dual-limb circuits

Dual-limb circuits require an active exhalation valve (p. 18) and an external flow sensor (p. 18).

Heated circuits for use with WILAméd AIRcon humidifier



WILAméd heated circuits are not available in North America.

Description	Quantity	Order number	Part number
Patient circuit, with dual heated limbs, with proximal pressure line, single-use, 22-mm ID (WILAméd)	10	989805661961	1127303
Patient circuit, with dual heated limbs, with proximal pressure line, single-use, 15-mm ID (WILAméd)	10	989805661971	1127304



Non-heated circuits

Description	Quantity	Order number	Part number
Patient circuit, dual-limb, non-heated, with proximal pressure line, 22-mm ID, single-use	10	989805661981	1127306
Patient circuit, dual-limb, non-heated, with proximal pressure line, 15-mm ID, single-use	10	989805661991	1127307



Heated circuits for use with Fisher & Paykel MR850 humidifier

Available with the purchase of a humidifier base only.


Description	Quantity	Order number	Part number
Patient circuit, with dual heated limbs, with pressure lines, single-use, adult (Fisher & Paykel RT111) (Not shown)	10	989805652821	1120221
Patient circuit, with dual heated limbs, with auto-feed chamber, with proximal pressure line, single-use, infant (Fisher & Paykel Evaqua RT385) (Not shown)	10	989805654581	1128391
Patient circuit, with dual heated limbs, with auto-feed chamber (flows > 4 L/min), with proximal pressure line, single-use, infant (Fisher & Paykel Evaqua RT265)	10	989805662901	1141548
Patient circuit, with dual heated limbs, with auto-feed chamber, with pressure lines, single-use, adult (Fisher & Paykel RT211)*	10	989805662891	1141427


* Exclusively available in the USA

Flex tracheal adapters

Adapters are 15 cm (6 in.) long. Adapters PN 1073902 and 1073904 connect to swivel and fixed, single-use exhalation ports.

Description	Quantity	Order number	Part number
Flex tracheal adapter, single-use, 22M/15F x 22F	10	989805633931	1073902
Flex tracheal adapter, single-use, 22M/15F x 22F	20	989805633951	1073904
Flex tracheal adapter, single-use, 22M/15F x 15M	10	989805633921	1073863
Flex tracheal adapter, single-use, 22M/15F x 15M	20	989805633941	1073903






Active exhalation valves

For single-limb use

Description	Quantity	Order number	Part number
Exhalation valve, active, single-limb (WILMed)	10	989805662061	1131661

For dual-limb use

Description	Quantity	Order number	Part number
Exhalation valve, active (AEV), dual-limb	10	989805662411	1139436
Exhalation valve, active (AEV), dual-limb	1	989805663541	1132110
Diaphragm, dual-limb active exhalation valve (AEV), dual-limb (Not shown)	10	989805662421	1139912



1139436
1132110

Bacterial/viral filters

The WILMed filter is not available in the USA.

Description	Quantity	Order number	Part number
Filter, bacterial/viral, single-use	10	989805609521	342077
Filter, bacterial/viral, single-patient use (WILMed)	10	989805662371	1135776



342077



1135776

External flow sensor

Description	Quantity	Order number	Part number
Flow sensor, with cable, external, Trilogy EV300, reusable, adult/pediatric	1	989805662071	1132106
Flow sensor, with cable, external, Trilogy EV300, reusable, pediatric/infant	1	989805662081	1132107
Flow sensor, external, Trilogy EV300, reusable, adult/pediatric	1	989805662351	1134711
Flow sensor, external, Trilogy EV300, reusable, pediatric/infant	1	989805662361	1134712
Flow sensor, external, Trilogy EV300, single-patient use, adult/pediatric	10	989805664651	1146497
Flow sensor, external, Trilogy EV300, single-patient use, pediatric/infant	10	989805664661	1146558
Cable, external flow sensor, Trilogy EV300	1	989805662321	1134592




1132106
1132107



1134711



1146497



1134712



1146558



1134592

FiO₂ sensor

Description	Order number	Part number
Sensor, FiO ₂	989805662341	1134668

External pulse oximeter

This highly accurate external pulse oximeter is for use with Nonin PureLight sensors.

Description	Quantity	Order number	Part number
Pulse oximeter, external (Nonin XPOD)	1	989805662051	1129640
Sensor, SpO ₂ , finger, reusable, adult (Nonin)	1	989805603171	936
Sensor, SpO ₂ , finger, reusable, pediatric (Nonin)	1	989805606971	936P
Sensor, SpO ₂ , cloth, disposable, adult (Nonin)	24	989805643671	1070692
Sensor, SpO ₂ , cloth, disposable, pediatric (Nonin)	24	989805643681	1070693
Sensor, SpO ₂ , cloth, disposable, infant (Nonin)	24	989805643691	1070694

Mainstream CO₂ sensor and accessories

Description	Quantity	Order number	Part number
Sensor, mainstream CO ₂	1	989803142651	M2501A
Airway adapter, mainstream CO ₂ , single-patient use, pediatric/adult (clear) (M2533)	10	989805618801	6063-00
Airway adapter, mainstream CO ₂ , single-patient use, neonatal (violet) (M2536A)	10	989805618811	6312-00
Airway adapter, mainstream CO ₂ , reusable, neonatal (red) (M2516A) (Not shown)	10	989805618921	7053-01
Airway adapter, mainstream CO ₂ , reusable, pediatric/adult (black) (M2513A) (Not shown)	1	989805618901	7007-01
Cable, adapter, mainstream CO ₂ sensor to USB	1	989805662031	1127404



7 Humidifiers

Fisher & Paykel MR850 humidifier

For details consult the Fisher & Paykel website, www.fphcare.com.

Humidifier bases

Model	Voltage	Languages	Power plug	Order number	Part number
AEK	230	English	United Kingdom (G)	989805656711	1134414
ADU	230	Danish	Denmark (K)	989805656721	1134416
AEU	230	English, French, Spanish, Greek, Czech, Croatian, Slovak, Romanian	Schuko (F)	989805656631	1134418
AFU	230	Dutch, French, German	Schuko (F)	989805656641	1134399
AGU	230	French, German, Italian	Schuko (F)	989805656691	1134404
ALU	230	Italian, Portuguese, Spanish	Schuko (F)	989805656241	1133848
ANU	230	Finnish, Swedish	Schuko (F)	989805656681	1134403
ARU	230	Bulgarian, Croatian, Czech, English, Hungarian, Polish, Slovak	Schuko (F)	989805656231	1133609
JHU	115	English, French, Spanish	USA (B)	989805634131	1084284



Humidifier chambers

Available with the purchase of a humidifier base only.

Description	Quantity	Order number	Part number
Humidifier chamber, Fisher & Paykel, single-use, adult (MR290)	10	989805642931	22104

Other humidification parts

Available with the purchase of a humidifier base only.

Description	Quantity	Order number	Part number
Adapter, heater wire, inspiratory-limb, MR850, for RT-series circuits (900MR806)	1	989805616981	1024934
Adapter, heater wire, dual (inspiratory/expiratory), MR850, for RT-series circuits (900MR805)	1	989805654511	1128070
Temperature probe, right-angle, for RT-series 1.83-m (6-ft) circuits (900MR861)	1	989805613101	1025039
Adapter, temperature probe, 22F x 22M/15F, single-use (900MR532)	1	989805642951	22161
Adapter, temperature probe, 22F x 22M, with 7.5-mm ID probe connector, single-use	50	453561517001	22164
Adapter, temperature probe, elbow, single-use (900MR131)	1	989805642961	22184



WILamed AIRcon humidifier

For details consult the WILamed website, <https://live.wilamed.de/?lang=en>

The WILamed humidifier is not available in North America.

Humidifier bases

Voltage	Languages	Power plug	Order number	Part number
230	English	United Kingdom (type G)	989805655781	1131790
230	French, German, Spanish	Continental Europe, (Schuko, type F)	989805654761	1129565
230	English	Israel (type H)	989805655801	1131792
230	Italian	Italy (type L)	989805655791	1131791
230	French, German, Italian	Switzerland (type J)	989805655811	1131793

Humidifier chambers

Description	Quantity	Order number	Part number
Humidifier chamber, WILamed AIRcon, auto-fill, single-use, adult	30	989805653111	1121825
Humidifier chamber, WILamed, single-use (500.380 C200AF Universal)	30	989805663041	1136694



8 Patient interfaces

Masks at a glance

NIV masks are classified according to intentional leak levels. Leak 1 masks have the lowest leak. Leak 1 masks require the use of a patient circuit with an exhalation port. SE masks are typically used for ventilation with a dual-limb circuit.

Model	Type	Patient population	Single-use	Reusable
PerforMax	Total face	7 year or older and > 20 kg	✓	✓
PerforMax pediatric (XXS and XS)	Total face	1 year or older and > 7 kg	✓	✓
AF531	Oro-nasal	7 years or older and > 20kg	✓	
AF541	Oro-nasal	7 years or older and > 20kg	✓	
PerformaTrak	Oro-nasal	> 30 kg	✓	
Contour Deluxe	Nasal	> 30 kg	✓	

Total face masks

Total face masks allow quick fitting and eliminate nasal bridge challenges by sealing around the perimeter of the face, where patients have less pressure sensitivity and smoother facial contours.

PerforMax total face mask, single-use

Interchangeable elbows available for the PerforMax easily adapt one mask to multiple ventilators. The masks shown below are intended for adult and pediatric patients (7 years or older and > 20 kg). For smaller pediatric patients, see the PerforMax pediatric ordering information on p. 32.

Order a global part unless a region-specific part is listed.

Elbow	Size	Quantity	Region	Order number	Part number
EE leak 2	S	10	Global	989805611281	1052535
EE leak 2	S	10	China	989805651481	CN1052535
EE leak 2	S	10	USA	989805611241	1052531
EE leak 2	L	10	Global	989805611291	1052536
EE leak 2	L	10	China	989805651491	CN1052536
EE leak 2	L	10	USA	989805611251	1052532
EE leak 2	XL	5	Global	989805634811	1080424
SE	S	10	Global	989805611301	1052547
SE	S	10	China	989805651501	CN1052547
SE	S	10	USA	989805611261	1052533
SE	L	10	Global	989805611311	1052548
SE	L	10	China	989805651511	CN1052548
SE	L	10	USA	989805611271	1052534
SE	XL	1	Global	989805634841	1080517
SE	XL	5	Global	989805634821	1080425



PerforMax single-use mask with EE leak 2 elbow



PerforMax single-use mask with SE elbow

PerforMax total face mask, reusable

Each mask includes one headgear. Some masks include multiple elbows. Order a global part unless a region-specific part is listed.

Elbow	Size	Quantity	Region	Order number	Part number
EE leak 2	L	1	USA	989805617091	1047427
EE leak 2	XL	1	USA	989805634761	1080419
EE leak 2 and SE	S	1	Global	989805611091	1044385
EE leak 2 and SE	L	1	Global	989805611101	1044386
EE leak 2 and SE	XL	1	Global	989805634801	1080423



PerforMax reusable mask with EE leak 2 elbow



PerforMax reusable mask with SE and EE leak 2 elbows

PerforMax total face mask set, reusable

Each mask set includes one mask and four headgear. Order a global part unless a region-specific part is listed.

Elbow	Size	Region	Order number	Part number
EE leak 2	S	Global	989805611381	1052555
EE leak 2	S	China	989805651551	CN1052555
EE leak 2	L	Global	989805611391	1052556
EE leak 2	L	China	989805651561	CN1052556
EE leak 2	XL	Global	989805634781	1080421
SE	S	Global	989805611401	1052568
SE	L	Global	989805611411	1052569
SE	XL	Global	989805634791	1080422
EE leak 2 and SE	S	Global	989805611361	1052553
EE leak 2 and SE	S	USA	989805611411	1049037
EE leak 2 and SE	L	Global	989805611371	1052554
EE leak 2 and SE	L	USA	989805611211	1049038
EE leak 2 and SE	XL	Global	989805634771	1080420



PerforMax mask set with EE leak 2 elbow



PerforMax mask set with SE elbow



PerforMax mask set with SE and EE leak 2 elbows

Oro-nasal masks

Oro-nasal masks are the standard of care for hospital NIV, and the Philips Respironics family includes a comprehensive selection of systems with features that address patient comfort, mask fitting, and ventilator compatibility.

AF541 oro-nasal mask, single-use

Designed to rest comfortably on the face, the AF541 mask features interchangeable under-the-nose and over-the-nose cushions. Alternating between the two cushion types can offload skin pressure points. Other features include a push-button forehead pad for easy mask adjustment and a nebulizer port for easy and efficient drug therapy administration. Interchangeable elbows can easily adapt one mask to multiple ventilators. The CapStrap headgear system simplifies mask application and routine patient care. The CleanClip system, available for the AF541, supports infection control; see p. 12.

The masks shown below are intended for adult and pediatric patients (7 years or older and > 20 kg). Order a global part unless a region-specific part is listed.

With over-the-nose cushion and CapStrap

The over-the-nose cushion, as in a traditional oro-nasal mask, sits on the bridge of the nose. Order a global part unless a region-specific part is listed.

Elbow	Size	Quantity	Region	Order number	Part number
EE leak 1 NE	S	10	Global	989805653931	1120925
EE leak 1 NE	S	10	USA	989805653851	1120915
EE leak 1 NE	M	10	Global	989805653941	1120926
EE leak 1 NE	M	10	USA	989805653861	1120916
EE leak 1 NE	L	10	Global	989805653951	1120927
EE leak 1 NE	L	10	USA	989805653871	1120917
EE leak 1 NE	XL	10	Global	989805653961	1120928
EE leak 1 NE	XL	10	USA	989805653881	1120918
EE leak 2 NE	S	10	Global	989805654011	1120935
EE leak 2 NE	M	10	Global	989805654021	1120936
EE leak 2 NE	L	10	Global	989805654031	1120937
EE leak 2 NE	XL	10	Global	989805654041	1120938
SE NE	S	10	Global	989805654091	1120945
SE NE	M	10	Global	989805654101	1120946
SE NE	L	10	Global	989805654111	1120947
SE NE	XL	10	Global	989805654121	1120948



AF541 mask with over-the-nose cushion, EE leak 1 NE elbow, and CapStrap



AF541 mask with over-the-nose cushion, EE leak 2 NE elbow, and CapStrap



AF541 mask with over-the-nose cushion, SE NE elbow, and CapStrap

With over-the-nose cushion and four-point headgear

The over-the-nose cushion, as in a traditional oro-nasal mask, sits on the bridge of the nose. Order a global part unless a region-specific part is listed.

Elbow	Size	Quantity	Region	Order number	Part number
EE leak 1 NE	S	10	Global	989805653891	1120920
EE leak 1 NE	S	10	USA	989805653811	1120910
EE leak 1 NE	M	10	Global	989805653901	1120921
EE leak 1 NE	M	10	USA	989805653821	1120911
EE leak 1 NE	L	10	Global	989805653911	1120922
EE leak 1 NE	L	10	USA	989805653831	1120912
EE leak 1 NE	XL	10	Global	989805653921	1120923
EE leak 1 NE	XL	10	USA	989805653841	1120913
EE leak 2 NE	S	10	Global	989805653971	1120930
EE leak 2 NE	M	10	Global	989805653981	1120931
EE leak 2 NE	L	10	Global	989805653991	1120932
EE leak 2 NE	XL	10	Global	989805654001	1120933
SE NE	S	10	Global	989805654051	1120940
SE NE	M	10	Global	989805654061	1120941
SE NE	L	10	Global	989805654071	1120942
SE NE	XL	10	Global	989805654081	1120943

With under-the-nose cushion and CapStrap

The under-the-nose cushion, which sits just below the nose, is offered in four sizes (XS, A, B, and C). Order a global part unless a region-specific part is listed.

Elbow	Size	Quantity	Region	Order number	Part number
EE leak 1 NE	XS	10	Global	989805656371	1133668
EE leak 1 NE	XS	10	USA	989805656291	1133660
EE leak 1 NE	A	10	Global	989805656381	1133669
EE leak 1 NE	A	10	USA	989805656301	1133661
EE leak 1 NE	B	10	Global	989805656391	1133670
EE leak 1 NE	B	10	USA	989805656311	1133662
EE leak 1 NE	C	10	Global	989805656401	1133671
EE leak 1 NE	C	10	USA	989805656321	1133663
EE leak 2 NE	XS	10	Global	989805656451	1133676
EE leak 2 NE	A	10	Global	989805656461	1133677
EE leak 2 NE	B	10	Global	989805656471	1133678
EE leak 2 NE	C	10	Global	989805656481	1133679
SE NE	XS	10	Global	989805656531	1133684
SE NE	A	10	Global	989805656541	1133685
SE NE	B	10	Global	989805656551	1133686
SE NE	C	10	Global	989805656561	1133687



AF541 mask with over-the-nose cushion, EE leak 1 NE elbow, and four-point headgear



AF541 mask with over-the-nose cushion, EE leak 2 NE elbow, and four-point headgear



AF541 mask with over-the-nose cushion, SE NE elbow, and four-point headgear



AF541 mask with under-the-nose cushion, EE leak 1 NE elbow, and CapStrap



AF541 mask with under-the-nose cushion, EE leak 2 NE elbow, and CapStrap



AF541 mask with under-the-nose cushion, SE NE elbow, and CapStrap

With under-the-nose cushion and four-point headgear

The under-the-nose cushion, which sits just below the nose, is offered in four sizes (XS, A, B, and C). Order a global part unless a region-specific part is listed.

Elbow	Size	Quantity	Region	Order number	Part number
EE leak 1 NE	XS	10	Global	989805656331	1133664
EE leak 1 NE	XS	10	USA	989805656251	1133626
EE leak 1 NE	A	10	Global	989805656341	1133665
EE leak 1 NE	A	10	USA	989805656261	1133627
EE leak 1 NE	B	10	Global	989805656351	1133666
EE leak 1 NE	B	10	USA	989805656271	1133658
EE leak 1 NE	C	10	Global	989805656361	1133667
EE leak 1 NE	C	10	USA	989805656281	1133659
EE leak 2 NE	XS	10	Global	989805656411	1133672
EE leak 2 NE	A	10	Global	989805656421	1133673
EE leak 2 NE	B	10	Global	989805656431	1133674
EE leak 2 NE	C	10	Global	989805656441	1133675
SE NE	XS	10	Global	989805656491	1133680
SE NE	A	10	Global	989805656501	1133681
SE NE	B	10	Global	989805656511	1133682
SE NE	C	10	Global	989805656521	1133683



AF541 mask with under-the-nose cushion, EE leak 1 NE elbow, and four-point headgear



AF541 mask with under-the-nose cushion, EE leak 2 NE elbow, and four-point headgear



AF541 mask with under-the-nose cushion, SE NE elbow, and four-point headgear

AF531 oro-nasal mask, single-use

Interchangeable elbows available for the AF531 easily adapt one mask to multiple ventilators. The CapStrap headgear system simplifies mask application and routine patient care. The CleanClip system, available for the AF531, supports infection control; see p. 12.

The masks shown below are intended for adult and pediatric patients (7 years or older and > 20 kg). Order a global part unless a region-specific part is listed.

With CapStrap headgear

Elbow	Size	Quantity	Region	Order number	Part number
EE leak 1	S	10	Global	989805628731	1072627
EE leak 1	S	10	China	989805652091	CN1072627
EE leak 1	M	1	China	989805652081	CN1061728
EE leak 1	M	10	Global	989805628741	1072628
EE leak 1	L	10	Global	989805628751	1072629
EE leak 1	L	10	China	989805652101	CN1072629
EE leak 2	S	10	Global	989805645621	1100390
EE leak 2	M	10	Global	989805645631	1100391
EE leak 2	L	10	Global	989805645641	1100392
SE	S	10	Global	989805628761	1072630
SE	M	10	Global	989805628771	1072631
SE	L	10	Global	989805628781	1072632
EE leak 1 NE	S	10	Global	989805634621	1072636
EE leak 1 NE	S	10	USA	989805635001	1072618
EE leak 1 NE	M	10	Global	989805634631	1072637
EE leak 1 NE	M	10	USA	989805635011	1072619
EE leak 1 NE	L	10	Global	989805634641	1072638
EE leak 1 NE	L	10	USA	989805635021	1072620
SE NE	S	10	Global	989805645681	1100396
SE NE	M	10	Global	989805645691	1100397
SE NE	L	10	Global	989805645701	1100438



AF531 mask with EE leak 1 elbow and CapStrap



AF531 mask with EE leak 2 elbow and CapStrap



AF531 mask with SE elbow and CapStrap

With four-point headgear

Elbow	Size	Quantity	Region	Order number	Part number
EE leak 1	S	10	Global	989805628671	1072621
EE leak 1	M	10	Global	989805628681	1072622
EE leak 1	M	1	China	989805652071	CN1061719
EE leak 1	L	10	Global	989805628691	1072623
EE leak 2	S	10	Global	989805645651	1100393
EE leak 2	S	10	China	989805649741	CN1100393
EE leak 2	M	10	Global	989805645661	1100394
EE leak 2	M	10	China	989805649751	CN1100394
EE leak 2	L	10	Global	989805645671	1100395
EE leak 2	L	10	China	989805649761	CN1100395
SE	S	10	Global	989805628701	1072624
SE	M	10	Global	989805628711	1072625
SE	M	10	China	989805649711	CN1072625
SE	L	10	Global	989805628721	1072626
EE leak 1 NE	S	10	Global	989805634591	1072633
EE leak 1 NE	S	10	USA	989805634971	1061733
EE leak 1 NE	M	10	Global	989805634601	1072634
EE leak 1 NE	M	10	USA	989805634981	1061734
EE leak 1 NE	L	10	Global	989805634611	1072635
EE leak 1 NE	L	10	USA	989805634991	1061735

AF531 mask with EE leak 1 elbow and four-point headgear



AF531 mask with EE leak 2 elbow and four-point headgear



AF531 mask with EE leak 1 NE elbow and four-point headgear



AF531 mask with SE elbow and four-point headgear



PerformaTrak oro-nasal mask, single-use

Dual-density foam forehead pads provide patient comfort and mask stability. Swivel clips simplify reapplication by eliminating the need to refit after removal. The PerformaTrak oro-nasal mask is intended for patients > 30 kg.

Order a global part unless a region-specific part is listed.

Elbow	Headgear	Size	Quantity	Region	Order number	Part number
EE leak 1	CapStrap	S	10	Global	989805610691	1019492
EE leak 1	CapStrap	S	10	USA	989805610631	1019478
EE leak 1	CapStrap	M	10	Global	989805610701	1019498
EE leak 1	CapStrap	M	10	USA	989805610641	1019479
EE leak 1	CapStrap	L	10	Global	989805610711	1019507
EE leak 1	CapStrap	L	10	USA	989805610651	1019480
EE leak 1	Four-point	S	10	Global	989805610481	1018589
EE leak 1	Four-point	S	10	USA	989805610511	1018592
EE leak 1	Four-point	M	10	Global	989805610491	1018590
EE leak 1	Four-point	M	10	USA	989805610521	1018593
EE leak 1	Four-point	L	10	Global	989805610501	1018591
EE leak 1	Four-point	L	10	USA	989805610531	1018594
SE	CapStrap	S	10	Global	989805610811	1019540
SE	CapStrap	S	10	USA	989805610751	1019511
SE	CapStrap	M	10	Global	989805610821	1019541
SE	CapStrap	M	10	USA	989805610761	1019525
SE	CapStrap	L	10	Global	989805610831	1019542
SE	CapStrap	L	10	USA	989805610771	1019533
SE	Four-point	S	10	Global	989805610581	1018605
SE	Four-point	S	10	USA	989805610551	1018596
SE	Four-point	M	10	Global	989805610591	1018606
SE	Four-point	M	10	USA	989805610561	1018597
SE	Four-point	L	10	Global	989805610541	1018595
SE	Four-point	L	10	USA	989805610571	1018598



PerformaTrak single-use oro-nasal mask with EE leak 1 elbow and CapStrap



PerformaTrak single-use oro-nasal mask with EE leak 1 elbow and four-point headgear



PerformaTrak single-use oro-nasal mask with SE elbow and CapStrap



PerformaTrak single-use oro-nasal mask with SE elbow and four-point headgear

Nasal masks

A selection of single-use nasal masks permit oral access and are particularly suitable for long-term noninvasive ventilation.

Contour Deluxe nasal mask, single-use

Dual-density foam forehead pads provide patient comfort and mask stability. Hook-and-loop grab tabs allow quick, easy removal and adjustment of headgear. The Contour Deluxe mask is intended for patients > 30 kg. Order a global part unless a region-specific part is listed.

Size	Quantity	Region	Order number	Part number
Petite	10	Global	989805610341	1016983
Petite	5	China	989805649691	CN1016980
Petite	5	USA	989805610311	1016980
S	10	Global	989805610301	1016698
S	10	USA	989805610261	1016694
S	5	Global	989805610291	1016697
S	5	China	989805649681	CN1016697
S	5	USA	989805610251	1016693
M/L	10	Global	989805610281	1016696
M/L	10	USA	989805610241	1016692
M/L	5	Global	989805610271	1016695
M/L	5	China	989805649671	CN1016695
M/L	5	USA	989805610231	1016691



Pediatric patient interfaces

PerforMax pediatric total face masks

Interchangeable elbows available for the PerforMax easily adapt one mask to multiple ventilators. The CleanClip system, available for the PerforMax, supports infection control; see p. 12.

Size XXS and XS interfaces are intended for pediatric patients 1 year or older and > 7 kg. Size S interfaces are intended for patients 7 years or older and > 20 kg. XXS and XS masks have child-friendly design features including colorful bonnets with a soft inner lining, color-matched clips, and a softer silicone material.

PerforMax pediatric total face mask with bonnet, single-use

Order a global part unless a region-specific part is listed.

Elbow	Size	Quantity	Region	Order number	Part number
EE leak 2	XXS	5	Global	989805645031	1083148
EE leak 2	XXS	5	China	989805651421	CN1083148
EE leak 2	XS	5	Global	989805645021	1083147
EE leak 2	XS	5	China	989805651441	CN1083147
SE	XXS	5	Global	989805645051	1083150
SE	XXS	5	China	989805651431	CN1083150
SE	XS	5	Global	989805645041	1083149
SE	XS	5	China	989805651461	CN1083149

PerforMax total face mask, single-use (size S)

Each mask includes a four-point headgear. Order a global part unless a region-specific part is listed.

Elbow	Size	Quantity	Region	Order number	Part number
EE leak 2	S	10	Global	989805611281	1052535
EE leak 2	S	10	USA	989805611241	1052531
SE	S	10	Global	989805611301	1052547
SE	S	10	USA	989805611261	1052533



PerforMax single-use pediatric mask with EE leak 2 elbow



PerforMax single-use pediatric mask with SE elbow



PerforMax single-use small mask with EE leak 2 elbow and four-point headgear



PerforMax single-use small mask with SE elbow and four-point headgear

PerforMax pediatric total face mask sets, reusable (sizes XXS and XS)

EE leak 2 mask sets include one mask, two each small pink and blue bonnets, and three additional elbow/tubing assemblies. SE mask sets include one mask and two each small pink and blue bonnets.

Elbow	Size	Order number	Part number
EE leak 2	XXS	989805644921	1083084
EE leak 2	XS	989805644911	1083083
SE	XXS	989805644941	1083086
SE	XS	989805644931	1083085



1083084
1083083

PerforMax total face mask, reusable (size S)

Each mask includes a four-point headgear. Order a global part unless a region-specific part is listed.

Elbow	Size	Region	Order number	Part number
EE leak 2	S	USA	989805617081	1047426
EE leak 2 and SE	S	Global	989805611091	1044385
EE leak 2 and SE	S	USA	989805611201	1049037



1083086
1083085

PerforMax total face mask set, reusable (size S)

Each mask set includes one mask and four each four-point headgear. Order a global part unless a region-specific part is listed.

Elbow	Size	Region	Order number	Part number
EE leak 2	S	Global	989805611381	1052555
SE	S	Global	989805611401	1052568
EE leak 2 and SE	S	Global	989805611361	1052553



1052555



1052568



1047426



1044385
1049037



1052553

AF531 oro-nasal mask, single-use (size S)

Interchangeable elbows available for the AF531 can easily adapt one mask to multiple ventilators. The CapStrap headgear system simplifies mask application and routine patient care. The CleanClip system, available for the AF531, supports infection control; see p. 12. This small version of the AF531 mask is intended for pediatric patients 7 years or older and > 20 kg. Order a global part unless a region-specific part is listed.

Elbow	Headgear	Size	Quantity	Region	Order number	Part number
EE leak 1	CapStrap	S	10	Global	989805628731	1072627
EE leak 1	CapStrap	S	10	China	989805652091	CN1072627
EE leak 1	Four-point	S	10	Global	989805628671	1072621
EE leak 2	CapStrap	S	10	Global	989805645621	1100390
EE leak 2	Four-point	S	10	Global	989805645651	1100393
EE leak 2	Four-point	S	10	China	989805649741	CN1100393
SE	CapStrap	S	10	Global	989805628761	1072630
SE	CapStrap	S	10	USA	989805642671	1094755
SE	Four-point	S	10	Global	989805628701	1072624
SE	Four-point	S	10	USA	989805642701	1094788
EE leak 1 NE	CapStrap	S	10	Global	989805634621	1072636
EE leak 1 NE	CapStrap	S	10	USA	989805635001	1072618
EE leak 1 NE	Four-point	S	10	Global	989805634591	1072633
EE leak 1 NE	Four-point	S	10	USA	989805634971	1061733
SE NE	CapStrap	S	10	Global	989805645681	1100396



AF531 mask with EE leak 1 elbow and CapStrap



AF531 mask with EE leak 1 elbow and four-point headgear



AF531 mask with SE elbow and CapStrap



AF531 mask with SE elbow and four-point headgear



AF531 mask with EE leak 2 elbow and CapStrap



AF531 mask with EE leak 2 elbow and four-point headgear



AF531 mask with EE leak 1 NE elbow and CapStrap



AF531 mask with SE NE elbow and CapStrap



AF531 mask with EE leak 1 NE elbow and four-point headgear

Pediatric nasal masks

PN841 pediatric nasal mask, single-patient use

Designed specifically for children, the reusable PN841 applies a child-friendly fabric pattern with a modified cushion curvature to provide a positive experience for the smallest patients. The leak correction dial on the PN841 enables clinicians to fix small leaks and readjust the fit of the cushion on the face. The PN841 is intended for patients > 10 kg (22 lb). Order a global part unless a region-specific part is listed.

Elbow	Size	Region	Quantity	Order number	Part number
SE	S	Global	5	989805655501	1130582
SE	S	USA Canada	5	989805655161	1126614
SE	M	Global	5	989805655511	1130583
SE	M	USA Canada	5	989805655171	1126615
SE	L	Global	5	989805655521	1130584
SE	L	USA Canada	5	989805655181	1126616



PN841 mask with SE elbow

9 Additional parts and accessories

Trilogy Remote Alarm

The Trilogy Remote Alarm is intended to call attention to an alarm condition with compatible Philips Respironics ventilators at a distance of up to 61 m (200 ft) from the device. An audible tone and a red, flashing indicator indicate that an alarm condition exists.

The Remote Alarm includes a 15-m (50-ft) RJ9 cable and a 9-V battery. The U.S version includes a power supply. If you are ordering the global version, also order a power supply.

Order a global part unless a region-specific part is listed.

Description	Order number	Part number
Remote Alarm, Trilogy, USA	989805653751	1118941
Remote Alarm, Trilogy, global	989805653761	1120632

Power supply for Remote Alarm

Order a global part unless a region-specific part is listed.

Description	Order number	Part number
Power supply, Trilogy Remote Alarm, USA	453561536261	1120633
Power supply, Trilogy Remote Alarm, global	453561537551	1120634
Power supply, Trilogy Remote Alarm, UK	453561537541	1127161

Remote Alarm stand

Description	Order number	Part number
Stand, Trilogy Remote Alarm (Not shown)	989805653741	34046

Remote Alarm cables

Description	Order number	Part number
Cable, Trilogy Remote Alarm (RJ9 to RJ9), 15-m (50-ft)	453561536241	1114933
Cable, Trilogy Remote Alarm (RJ9 to RJ9), 61-m (200-ft)	453561537551	1114934



1118941
1120632



1120633
1120634
1127161



1114933
1114934


Operator replacement parts

Description	Quantity	Order number	Part number
Battery, Trilogy EV300, detachable, global	1	989805661831	1127889
Battery, Trilogy EV300, detachable, India	1	989805662401	1136881
Filter, air inlet, foam, Trilogy EV300	10	989805662311	1134591
Filter, particulate, Trilogy EV300	10	989805662301	1134430
Clamp, power cord, Trilogy EV300 (Not shown)	1	989805662271	1134347
Quick-connect insert, O ₂ , Trilogy	2	989805647691	1040390



Nurse call cables

Description	Order number	Part number
Cable, nurse call, RJ9 to 1/4-in. mono jack, normally open (alarm state = closed), 3.65-m (12-ft)	989805635491	1080247
Cable, nurse call, RJ9 to 1/4-in. mono jack, normally closed (alarm state = open), 3.65-m (12-ft)	989805635501	1080248
Cable, nurse call, RJ9 to open end, 3-m (9.8-ft)	989805635511	1080249



Power cables

Description	Order number	Part number
Cable, 12/24-V battery, with terminals	989805662001	1127401
Cable, 12/24-V battery, with car adapter	989805662011	1127402



Patient monitor interface cables

Not shown.

Description	Order number	Part number
Cable, USB to DB9, isolated. Interfaces the EV300 with a hospital monitor such as the Alice 6 PSG System, or with a computer that is running Omnilab Direct.	989805662041	1127405
Cable, USB to RJ45, isolated. Interfaces the EV300 with a hospital monitor such as the Philips IntelliVue.	989805662021	1127403

Carrying case

The carrying case may be used while the ventilator is operating.

Description	Order number	Part number
Carrying case, Trilogy EV300, in-use	989805662261	1133930



1133930

Service manual

Description	Order number	Part number
Service manual, Trilogy EV300 (English) (Not shown)	989805661931	1144577

Test lung

A test lung is required for ventilator testing. 1-liter capacity. Latex free.

Description	Order number	Part number
Test lung, hard-sided (low-compliance)	989805611871	1021671



1021671

Part number index

Numbers

936 19	1048899 13	1080423 24	1121825 21	1133674 27
936P 19	1049037 24, 33	1080424 23	1126614 35	1133675 27
6063-00 19	1049038 24	1080425 23	1126615 35	1133676 26
6312-00 19	1052531 23, 32	1080517 23	1126616 35	1133677 26
7007-01 19	1052532 23	1083083 33	1127161 36	1133678 26
7053-01 19	1052533 23, 32	1083084 33	1127303 16	1133679 26
22104 20	1052534 23	1083085 33	1127304 16	1133680 27
22161 20	1052535 23, 32	1083086 33	1127306 16	1133681 27
22164 20	1052536 23	1083147 32	1127307 16	1133682 27
22184 20	1052547 23, 32	1083148 32	1127401 37	1133683 27
34046 36	1052548 23	1083149 32	1127402 37	1133684 26
311107 12	1052553 24, 33	1083150 32	1127403 37	1133685 26
332497 12	1052554 24	1084284 20	1127404 19	1133686 26
342077 18	1052555 24, 33	1094755 34	1127405 37	1133687 26
1001659 9	1052556 24	1094788 34	1127889 37	1133848 20
1001663 9	1052568 24, 33	1100390 28, 34	1128070 20	1133930 38
1001664 7, 9	1052569 24	1100391 28	1128391 17	1134347 6, 37
1002229 13	1052863 13	1100392 28	1129565 21	1134399 20
1002230 13	1054832 7, 9	1100393 29, 34	1129640 19	1134403 20
1002497 13	1054833 7, 9	1100394 29	1130582 35	1134404 20
1002506 13	1054834 7, 9	1100395 29	1130583 35	1134414 20
1002556 9	1054835 9	1100396 28, 34	1130584 35	1134416 20
1003661 12	1054836 7, 9	1100397 28	1131463 11	1134418 20
1009753 13	1061733 29, 34	1100438 28	1131661 18	1134430 6, 37
1016691 31	1061734 29	1109602 11	1131790 21	1134591 37
1016692 31	1061735 29	1109881 11	1131791 21	1134592 18
1016693 31	1062978 7, 9	1113371 11	1131792 21	1134668 19
1016694 31	1070692 19	1113392 11	1131793 21	1134711 18
1016695 31	1070693 19	1114933 36	1132041 11	1134712 18
1016696 31	1070694 19	1114934 36	1132106 18	1135776 18
1016697 31	1071789 9	1118941 36	1132107 18	1136694 21
1016698 31	1071790 7, 9	1120221 17	1132110 18	1136845 8
1016980 31	1071791 7, 9	1120632 36	1132339 14	1136881 37
1016983 31	1071792 9	1120633 36	1132340 6, 14	1139436 18
1018589 30	1072618 28, 34	1120634 36	1132341 14	1139909 14
1018590 30	1072619 28	1120910 26	1132342 14	1139912 18
1018591 30	1072620 28	1120911 26	1132343 15	1141427 17
1018592 30	1072621 29, 34	1120912 26	1132344 15	1141548 17
1018593 30	1072622 29	1120913 26	1132345 15	1141764 11
1018594 30	1072623 29	1120915 25	1132346 15	1141765 11
1018595 30	1072624 29, 34	1120916 25	1133284 7, 8	1141766 11
1018596 30	1072625 29	1120917 25	1133286 8	1143216 10
1018597 30	1072626 29	1120918 25	1133287 7, 8	1143217 7, 10
1018598 30	1072627 28, 34	1120920 26	1133468 7, 8	1143278 7, 10
1018605 30	1072628 28	1120921 26	1133469 7, 8	1143279 7, 10
1018606 30	1072629 28	1120922 26	1133470 8	1143280 7, 10
1019478 30	1072630 28, 34	1120923 26	1133471 7, 8	1143281 10
1019479 30	1072631 28	1120925 25	1133472 7, 8	1143282 7, 10
1019480 30	1072632 28	1120926 25	1133609 20	1143284 10
1019492 30	1072633 29, 34	1120927 25	1133626 27	1143285 7, 10
1019498 30	1072634 29	1120928 25	1133627 27	1143286 10
1019507 30	1072635 29	1120930 26	1133658 27	1143287 7, 10
1019511 30	1072636 28, 34	1120931 26	1133659 27	1143288 7, 10
1019525 30	1072637 28	1120932 26	1133660 26	1143289 7, 10
1019533 30	1072638 28	1120933 26	1133661 26	1143290 10
1019540 30	1073863 17	1120935 25	1133662 26	1143291 10
1019541 30	1073902 17	1120936 25	1133663 26	1143292 10
1019542 30	1073903 17	1120937 25	1133664 27	1143293 7, 10
1021671 38	1073904 17	1120938 25	1133665 27	1143295 7, 10
1024934 20	1076430 12	1120940 26	1133666 27	1143296 10
1025039 20	1080247 37	1120941 26	1133667 27	1143297 10
1040390 37	1080248 37	1120942 26	1133668 26	1143298 10
1044385 24, 33	1080249 37	1120943 26	1133669 26	1143299 7, 10
1044386 24	1080419 24	1120945 25	1133670 26	1143300 7, 10
1047426 33	1080420 24	1120946 25	1133671 26	1143301 10
1047427 24	1080421 24	1120947 25	1133672 27	1143302 7, 10
	1080422 24	1120948 25	1133673 27	1143303 7, 10

1143304 7, 10	989805611391 24	989805642671 34	989805653961 25	989805660781 11
1143305 7, 10	989805611401 24, 33	989805642701 34	989805653971 26	989805660791 11
1143306 7, 10	989805611411 24	989805642931 20	989805653981 26	989805661831 37
1143307 10	989805611511 13	989805642951 20	989805653991 26	989805661841 5
1143308 10	989805611871 38	989805642961 20	989805654001 26	989805661851 5
1143309 7, 10	989805612421 9	989805643161 14	989805654011 25	989805661861 5
1143310 7, 10	989805612551 7, 9	989805643671 19	989805654021 25	989805661871 5
1143311 7, 10	989805612691 9	989805643681 19	989805654031 25	989805661881 5
1143312 10	989805612701 9	989805643691 19	989805654041 25	989805661891 5
1144577 38	989805613101 20	989805644911 33	989805654051 26	989805661911 5
1146497 18	989805613131 13	989805644921 33	989805654061 26	989805661931 38
1146558 18	989805613291 13	989805644931 33	989805654071 26	989805661961 16
453561506591 13	989805613311 13	989805644941 33	989805654081 26	989805661971 16
453561507401 12	989805613541 13	989805645021 32	989805654091 25	989805661981 16
453561517001 20	989805613851 7, 9	989805645031 32	989805654101 25	989805661991 16
453561536241 36	989805615961 7, 9	989805645041 32	989805654111 25	989805662001 37
453561536261 36	989805615971 7, 9	989805645051 32	989805654121 25	989805662011 37
453561537541 36	989805616201 7, 9	989805645621 28, 34	989805654511 20	989805662021 37
453561537551 36	989805616871 9	989805645631 28	989805654581 17	989805662031 19
989803142651 19	989805616881 7, 9	989805645641 28	989805654761 21	989805662041 37
989805603171 19	989805616981 20	989805645651 29, 34	989805655161 35	989805662051 19
989805605731 13	989805617081 33	989805645661 29	989805655171 35	989805662061 18
989805606971 19	989805617091 24	989805645671 29	989805655181 35	989805662071 18
989805609521 18	989805617961 12	989805645681 28, 34	989805655501 35	989805662081 18
989805610231 31	989805618801 19	989805645691 28	989805655511 35	989805662091 14
989805610241 31	989805618811 19	989805645701 28	989805655521 35	989805662101 6, 14
989805610251 31	989805618901 19	989805647691 37	989805655781 21	989805662111 14
989805610261 31	989805618921 19	989805648241 11	989805655791 21	989805662121 14
989805610271 31	989805621681 9	989805648741 11	989805655801 21	989805662131 15
989805610281 31	989805621691 7, 9	989805648891 12	989805655811 21	989805662141 15
989805610291 31	989805621701 7, 9	989805649671 31	989805656061 11	989805662151 15
989805610301 31	989805621711 9	989805649681 31	989805656081 11	989805662161 15
989805610311 31	989805627181 12	989805649691 31	989805656231 20	989805662171 7, 8
989805610341 31	989805628671 29, 34	989805649711 29	989805656241 20	989805662191 8
989805610481 30	989805628681 29	989805649741 29, 34	989805656251 27	989805662201 7, 8
989805610491 30	989805628691 29	989805649751 29	989805656261 27	989805662211 7, 8
989805610501 30	989805628701 29, 34	989805649761 29	989805656271 27	989805662221 7, 8
989805610511 30	989805628711 29	989805650181 11	989805656281 27	989805662231 8
989805610521 30	989805628721 29	989805650191 11	989805656291 26	989805662241 7, 8
989805610531 30	989805628731 28, 34	989805651421 32	989805656301 26	989805662251 7, 8
989805610541 30	989805628741 28	989805651431 32	989805656311 26	989805662261 38
989805610551 30	989805628751 28	989805651441 32	989805656321 26	989805662271 6, 37
989805610561 30	989805628761 28, 34	989805651461 32	989805656331 27	989805662301 6, 37
989805610571 30	989805628771 28	989805651481 23	989805656341 27	989805662311 37
989805610581 30	989805628781 28	989805651491 23	989805656351 27	989805662321 18
989805610591 30	989805633921 17	989805651501 23	989805656361 27	989805662341 19
989805610631 30	989805633931 17	989805651511 23	989805656371 26	989805662351 18
989805610641 30	989805633941 17	989805651551 24	989805656381 26	989805662361 18
989805610651 30	989805633951 17	989805651561 24	989805656391 26	989805662371 18
989805610691 30	989805634131 20	989805652071 29	989805656401 26	989805662381 8
989805610701 30	989805634591 29, 34	989805652081 28	989805656411 27	989805662401 37
989805610711 30	989805634601 29	989805652091 28, 34	989805656421 27	989805662411 18
989805610751 30	989805634611 29	989805652101 28	989805656431 27	989805662421 18
989805610761 30	989805634621 28, 34	989805652821 17	989805656441 27	989805662431 10
989805610771 30	989805634631 28	989805653111 21	989805656451 26	989805662441 7, 10
989805610811 30	989805634641 28	989805653741 36	989805656461 26	989805662451 7, 10
989805610821 30	989805634761 24	989805653751 36	989805656471 26	989805662461 7, 10
989805610831 30	989805634771 24	989805653761 36	989805656481 26	989805662471 7, 10
989805611091 24, 33	989805634781 24	989805653811 26	989805656491 27	989805662481 10
989805611101 24	989805634791 24	989805653821 26	989805656501 27	989805662491 7, 10
989805611201 33	989805634801 24	989805653831 26	989805656511 27	989805662511 10
989805611211 24	989805634811 23	989805653841 26	989805656521 27	989805662521 7, 10
989805611241 23, 32	989805634821 23	989805653851 25	989805656531 26	989805662531 10
989805611251 23	989805634841 23	989805653861 25	989805656541 26	989805662541 7, 10
989805611261 23, 32	989805634971 29, 34	989805653871 25	989805656551 26	989805662551 7, 10
989805611271 23	989805634981 29	989805653881 25	989805656561 26	989805662561 7, 10
989805611281 23, 32	989805634991 29	989805653891 26	989805656631 20	989805662571 10
989805611291 23	989805635001 28, 34	989805653901 26	989805656641 20	989805662581 10
989805611301 23, 32	989805635011 28	989805653911 26	989805656681 20	989805662591 10
989805611311 23	989805635021 28	989805653921 26	989805656691 20	989805662601 7, 10
989805611361 24, 33	989805635491 37	989805653931 25	989805656711 20	989805662621 7, 10
989805611371 24	989805635501 37	989805653941 25	989805656721 20	989805662631 10
989805611381 24, 33	989805635511 37	989805653951 25	989805660771 11	989805662641 10

989805662651 10	GB2200X15B 5
989805662661 7, 10	
989805662671 7, 10	IA2200X15B 5
989805662681 10	IN2200X15B 5
989805662691 7, 10	IT2200X21B 5
989805662701 7, 10	
989805662711 7, 10	KR2200X15B 5
989805662721 7, 10	
989805662731 7, 10	M2501A 19
989805662741 10	
989805662751 10	ND2200X15B 5
989805662761 7, 10	
989805662771 7, 10	RU2200X15B 5
989805662781 7, 10	
989805662791 10	TR2200X15B 5
989805662891 17	
989805662901 17	
989805663001 5	
989805663011 5	
989805663041 21	
989805663381 5	
989805663391 5	
989805663401 5	
989805663411 5	
989805663421 5	
989805663431 5	
989805663441 5	
989805663451 5	
989805663461 5	
989805663541 18	
989805664651 18	
989805664661 18	
989805667221 5	
BL2200X15B 5	
BR2200X18B 5	
CA2200X12B 5	
CN2200X17B 5	
CN105235 23	
CN1016695 31	
CN1016697 31	
CN1016980 31	
CN1052535 23	
CN1052536 23	
CN1052547 23	
CN1052548 23	
CN1052555 24	
CN1052556 24	
CN1061719 29	
CN1061728 28	
CN1072625 29	
CN1072627 28, 34	
CN1072629 28	
CN1083147 32	
CN1083148 32	
CN1083149 32	
CN1083150 32	
CN1100393 29, 34	
CN1100394 29	
CN1100395 29	
DE2200X13B 5	
DS2200X11B 5	
EE2200X15B 5	
ES2200X15B 5	
EU2200X15B 5	
EU2200X19 5	
FR2200X14B 5	
FX2200X15B 5	

Alphabetic index

A

active exhalation valves 18
AF531 oro-nasal mask
 adult 28
 pediatric 34
AF541 oro-nasal mask 25
airway adapters, CO₂ sensor 19
arm, support. *See* support arm

B

battery, detachable external 37
brackets, mounting
 CleanClip mask storage system 12
 support arm 12
breathing circuits. *See* patient circuits

C

cables 37
carrying case 38
cart. *See* stand, roll
circuits, patient. *See* patient circuits
CleanClip mask storage system 12
CO₂ sensor and accessories 19
Contour Deluxe nasal mask 31
cylinder holder, O₂ 11

E

exhalation port, swivel 14
exhalation valves, active 18

F

face masks. *See* total face masks, oro-nasal masks, or name of specific mask
filter
 air inlet 37
 bacterial/viral 18
 particulate 37
FiO₂ sensor 19
Fisher & Paykel MR850 humidifier and accessories 20
fittings, O₂ 9
flex arm. *See* support arm
flow sensor 18
full face masks. *See* oro-nasal masks or name of specific mask

H

hanger, tubing 12
hoses, O₂ 9
humidifiers 20
 Fisher & Paykel MR850 20
 WILMed AIRcon 21

I

idle mask storage bracket 12

M

manuals, technical 10, 38
masks. *See* patient interfaces or name of specific mask
mounting brackets. *See* brackets, mounting
MR850 humidifier. *See* Fisher & Paykel MR850 humidifier and accessories

N

nasal masks. *See also* name of specific mask
 adult 32
 pediatric 35

O

O₂ cylinder holder. *See* cylinder holder, O₂
O₂ fittings 9
O₂ hoses 9
O₂ manifold/transport kit 11
O₂ sensor 19
operator replacement parts 37
oro-nasal masks 25. *See also* name of specific mask

P

patient circuits 14–19
patient interfaces 22–32
 AF531 oro-nasal mask 28, 34
 AF541 oro-nasal mask 25
 Contour Deluxe nasal mask 31
 nasal masks
 adult 31
 pediatric 35
 oro-nasal masks
 adult 25
 pediatric 32
 PerformaTrak oro-nasal mask
 single-use 30
 PerforMax total face mask
 reusable, adult 24
 reusable, pediatric 33
 single-use, adult 23
 single-use, pediatric 32
 PN841 pediatric nasal mask 35
 total face masks
 adult 23
 pediatric 32
 PerformaTrak nasal mask 31
 PerformaTrak oro-nasal mask 30
 single-use 30
 PerforMax total face mask
 reusable
 adult 24
 pediatric 33
 single-use
 adult 23
 pediatric 32
 PN841 pediatric nasal mask 35
power cords 8
pulse oximeter and accessories 19

R

Remote Alarm 36
roll stand 11

S

sensor
 CO₂ 19
 FiO₂ 19
 flow 18
service manual 38
small product roll stand (SPRS-2) conversion kit 11

support arm 12
swivel exhalation port 14

T

technical manuals 10
test lung 38
total face masks
 adult 23
 pediatric 32
transport kit, O₂ 11
Trilogy Remote Alarm 36
trolley. *See* stand, roll
tubing circuits. *See* patient circuits
tubing hanger 12

U

universal stand. *See*
user manuals 10

V

valves, exhalation, active 18
ventilator circuits. *See* patient circuits

W

water bag holders 13
WILMed AIRcon humidifier 21

