



PHILIPS

Image guided therapy

Zenition 70 series

Touch screen module

Zenition 70 with
Touch screen module
Unlimited potential at your fingertips

More imaging control at table side

The Touch screen module (TSM) is designed to deliver value by enhancing the clinical capabilities of Zenition 70. As an addition to the Zenition 70 platform, the TSM can contribute to increased surgeon table-side control during procedures. The optional touch screen boasts the same functionality as the C-arm-mounted screen and displays the same content as the Mobile View Station.



Exceptional imaging clarity and flexibility

1. Flexibility

Control imaging parameters from within the sterile field

2. Positioning

Adjust position through multiple angles with the swing arm

3. Ease of use

With a large 12.1-inch display, the TSM provides similar ease of use to a stand UI

4. Performance

Instantaneous touch response (from TSM to stand user interface)



Hassle-free communications, greater control

The TSM lets you control C-arm functions directly on the touchscreen – allowing you to select, drag and zoom items with ease. This means you can view images instantly, giving you greater image control during procedures. Fast access to images on the TSM paves the way for smooth, straightforward and efficient communication.

Touch screen module

Key specifications

Specifications TSM:

Size	12.1 inch LCD
Position	Rotate and tilt
Display matrix	1280 x 800
Weight	2.05 kg
Dimensions	Width 310.75 mm
	Depth 53.84 mm
	Height 212.75 mm
Cable length	5 meters
Environmental conditions classification according to IEC60529	IP23

The TSM solution comprises:

Swing arm, breakaway connector and Touch screen monitor

Please contact your Philips representative for local availability of the Touch screen module.

